



BIOLABO
 www.biolabo.fr
MANUFACTURER:
BIOLABO SAS,
 Les Hautes Rives
 02160, Maizy, France

HAEMOGLOBIN

Colorimetric Method (Cyanmethemoglobin)

Reagent for quantitative determination of haemoglobin (Hb)
 in whole blood

REF 3502200 R1 2 x 200 mL



IVD IN VITRO DIAGNOSTIC USE

TECHNICAL SUPPORT AND ORDERS

Tel : (33) 03 23 25 15 50

Fax: (33) 03 23 256 256

CLINICAL SIGNIFICANCE (1)

A sufficient concentration of haemoglobin in blood is essential for adequate transport of O₂ and CO₂ between lungs and other tissues. Blood haemoglobin concentration may be diminished as a consequence of hemorrhage or hemolysis or as a result of impaired blood formation in bone marrow. Conversely, blood haemoglobin concentration may be increased when gas exchange through the lungs is impaired or in various other disorders. Measurement of the blood haemoglobin concentration is important as an initial step in the detection of anemia (diminished haemoglobin concentration) or erythrocytosis (increased red blood cells count and haemoglobin concentration).

PRINCIPLE (4) (5)

Method recognised as reference method by ICSH (International Committee of Standardisation in Hematology).

Fe²⁺ of haemoglobin is oxidised to the Fe³⁺ of methaemoglobin by ferricyanide, and the methaemoglobin is converted into stable cyanmethaemoglobin by addition of potassium cyanide (KCN).

The absorbance of cyanmethaemoglobin, directly proportional to the haemoglobin concentration, is measured at 546 nm (520-560).

REAGENT COMPOSITION

Vial R1 WORKING REAGENT

Phosphate Buffer	1	mmol/L
Potassium cyanide	0,75	mmol/L
Potassium ferricyanide	0,6	mmol/L
Detergent	0.1	g/L
Preservative	< 0.1	%

SAFETY CAUTIONS

BIOLABO reagents are designated for professional, in vitro diagnostic use.

- Verify the integrity of the contents before use.
- Use adequate protections (overall, gloves, glasses).
- Do not pipette by mouth.
- In case of contact with skin and eyes, thoroughly wash affected areas with plenty of water and seek medical advice.
- Material Safety Data Sheet is available upon request.
- Waste disposal: Respect legislation in force in the country.

All specimens should be handled as potentially infectious, in accordance with good laboratory practices using appropriate precautions. Respect legislation in force in the country.

REAGENT PREPARATION

Reagent is ready for use.

STABILITY AND STORAGE

Store away from light, well cap in the original vial at 18-25°C.

- Reagents are stable until expiry date stated on the label of the kit when stored and used as described in the insert and free from contamination.
- Discard any reagent if cloudy or if reagent blank at 546 nm is > 0.010.

SPECIMEN COLLECTION AND HANDLING (2)

Whole blood (EDTA).

Fetal blood: collect by percutaneous umbilical blood sampling.

Swirl gently to homogenise before assay.

Haemoglobin is stable in specimen for:

- 48 h at 2-8°C.
- 24 h at room temperature (< 25°C).

INTERFERENCES (2) (3)

Lipemia or leucocitar concentration > 25.10⁹/L involve overestimated results. Overestimation have been detected in the presence of HbC or HbS, in serious liver disorders or in globulin precipitation (ex: multiple myelome or Waldenström macroglobulinemia).

For a more comprehensive review of factors affecting this assay refer to the publication of Young D.S.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Basic medical analysis laboratory equipment.
2. Normal and pathological control blood.

CALIBRATION (6)

Use the calibration factor indicated in § **CALCULATION** or a calibrator (cyanmethemoglobin) assayed with the same method.

The calibration frequency depends on proper instrument functions and on preservation of the reagent.

It is recommended to calibrate in the following cases :

1. When changing vial of reagent.
2. After maintenance operations on the instrument.
3. If control values are out of range, even after using a new vial of fresh blood



Manufacturer



Use by



In vitro diagnostic



Temperature limitation



Catalogue number



See insert



Batch number



Store away from light



sufficient for



dilute with

QUALITY CONTROL

- Assayed control blood referring to the same method.
- External quality control program.

It is recommended to control in the following cases:

- At least once a run.
- At least once within 24 hours.
- When changing vial of reagent.
- After maintenance operations on the instrument.

If control is out of range, apply following actions:

1. Repeat the test with the same control.
2. If control is still out of range, prepare a fresh control and repeat the test.
3. If control is still out of range, verify analysis parameters: Wavelength, specimen/reagent ratio, calibration factor.
4. If control is still out of range, use a new vial of reagent and re-assay.
5. If control is still out of range, please contact BIOLABO technical support or your local Agent.

EXPECTED VALUES (2)

In fetal blood	g/dL	g/L	mmol/L
18-20 weeks	11.5 ± 0.78	115 ± 7.8	7.13 ± 0.48
21-22 weeks	12.3 ± 0.89	123 ± 8.9	7.63 ± 0.55
23-25 weeks	12.4 ± 0.77	124 ± 7.7	7.69 ± 0.48
26-30 weeks	13.4 ± 1.17	134 ± 12	8.31 ± 0.75

In cord blood	g/dL	g/L	mmol/L
	13.5-20.5	135-205	8.37-12.7

In total blood	g/dL	g/L	mmol/L
0.5 months	13.4-19.8	134-198	8.31-12.28
1 months	10.7-17.1	107-171	6.63-10.6
2 months	9.4-13.0	94-130	5.83-8.06
4 months	10.3-14.1	103-141	6.39-8.74
6 months	11.1-14.1	111-141	6.88-8.74
9 months	11.4-14.0	114-140	7.07-8.68
12 months	11.3-14.1	113-141	7.01-8.74
1-2 years	11.0-14.0	110-140	6.82-8.68
2-5 years	11.0-14.0	110-140	6.82-8.68
5-9 years	11.5-14.5	115-145	7.13-8.99
9-12 years	12.0-15.0	120-150	7.44-9.3
12-14 years	M 12.0-16.0	120-160	7.44-9.92
	F 11.5-15.0	115-150	7.13-9.3
15-17 years	M 11.7-16.6	117-166	7.25-10.29
	F 11.7-15.3	117-153	7.25-9.49
18-44 years	M 13.2-17.3	132-173	8.18-10.73
	F 11.7-15.5	117-155	7.25-9.61
45-64 years	M 13.1-17.2	131-172	8.12-10.66
	F 11.7-16.0	117-160	7.25-9.92
65-74 years	M 12.6-17.4	126-174	7.81-10.79
	F 11.7-16.1	117-161	7.25-9.98

It is recommended that each laboratory establish its own normal ranges for the population that it serves.

PERFORMANCES

Within run N = 20	Low level	High level	Between run N = 20	Low level	High level
Mean g/dL	6.7	18.9	Mean g/dL	6.3	17.1
S.D. g/dL	0.05	0.1	S.D. g/dL	0.29	0.42
C.V. %	0.7	0.6	C.V. %	4.6	2.5

Detection limit: approximately 0.3 g/dL

Sensitivity for 10 g/dL: approximately 0.272 Abs at 546 nm.

Comparison study with commercially available reagent:

$$y = 0.9999x + 0.08$$

$$r = 0.9962$$

LINEARITY

The assay is linear up to 250 g/L, 25 g/dL, 15.5 mmol/L (Hb/4).

MANUAL PROCEDURE

Pipette into test tubes.	Blank	Assay
Reagent	5 mL	5 mL
Demineralised water	20 µL	
Homogenized blood		20 µL

It is recommended to use a positive moved pipette to dispense blood. Rinse pipette several times into the reagent. Mix well and incubate at least for 3 minutes at room temperature. Read absorbance at 546 nm (520-560) against reagent blank.

Away from light, reaction is stable at least for 1 hour.

Note: Specific procedures are available upon request for automated instruments. Please contact BIOLABO technical support.

CALCULATION

	λ = 530 nm	λ = 546 nm	λ = 550 nm
Hb (g/L)	Abs x 386.1	Abs x 367.7	Abs x 376.2
Hb (g/dL)	Abs x 38.61	Abs x 36.77	Abs x 37.62
Hb mmol/L (Hb/4)	Abs x 23.96	Abs x 22.82	Abs x 23.34

This factors were designated as a guide only and may slightly vary. It is recommended to verify with control blood.

REFERENCES

- (1) TIETZ N.W. *Text book of clinical chemistry*, 3rd Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p. 1657-1688.
- (2) *Clinical Guide to Laboratory Test*, 3rd Ed., N.W. TIETZ (1995) p. 312-314.
- (3) YOUNG D.S., *Effect of Drugs on Clinical laboratory Tests*, 4th Ed. (1995) p. 3-325 to 3-330
- (4) DRABKIN, D.L., and AUSTIN, J.H., *J Biol. Chem.*, (1935), 112, p.51
- (5) VAN KAMPEN, E.J. and ZIJLSTRA W.G., *Determination of haemoglobin and its derivatives advances in clinical chemistry* (1965), 8, 141-187
- (6) VAN KAMPEN, E.J. and ZIJLSTRA W.G., *International Committee for standardization in haematology, British journal of haematology* (1967), 13 [Suppl] 71