



BIOLABO
www.biolabo.fr

MANUFACTURER:
BIOLABO SAS,

Les Hautes Rives
02160, Maizy, France

Inorganic PHOSPHORUS

U.V.Method

Reagent for quantitative determination of inorganic phosphate
in human plasma and serum, or urines

REF 80015 R1 2 x 125 mL R2 1 x 5 mL

TECHNICAL SUPPORT AND ORDERS

Tel: (33) 03 23 25 15 50

Fax: (33) 03 23 256 256



IVD IN VITRO DIAGNOSTIC USE

CLINICAL SIGNIFICANCE (1) (2)

An adult human body contains approximately 600 g of phosphates expressed as phosphorus, of which about 85% is bound to calcium in bones and the rest principally in other tissue cells. Most phosphate present in cells is organic and incorporated into phospholipids, nucleic acids, and high energy compounds. Serum/plasma contains approximately 1% of total phosphate as inorganic phosphate, the fraction measured in routine biochemical analysis.

An elevation of phosphorus in serum/plasma is often associated with bone diseases, renal failures, hypoparathyroidism, hypervitaminosis D...

Decreased serum/plasma phosphorus concentrations are found in case of osteomalacia, vitamin D deficiency, primary hyperparathyroidism...

PRINCIPLE (4) (5)

Method without deproteinisation described by Daly and al. modified by Gamst O.K. and Try K.

In an acid medium, phosphate ions form a phosphomolybdic complex with the ammonium molybdate. The absorbance measured at 340 nm is proportional to the concentration of phosphate ions in the specimen.

REAGENTS COMPOSITION

Vial R1

MOLYBDATE REAGENT

Ammonium Molybdate 0.63 mmol/L
Sulfuric acid 210 mmol/L
Surfactant

Xi, R36/38: Irritating to eyes and skin

S36/37/39: Wear suitable protective clothing, gloves and eye/face protection

S26-S28: In case of contact of eyes, rinse immediately with plenty of water and seek medical advice. In case of contact with skin, rinse immediately with plenty of water.

Vial R2

STANDARD

Phosphorus 5 mg/dL (1.61 mmol/L)

SAFETY CAUTIONS

BIOLABO reagents are designated for professional, in vitro diagnostic use.

- Verify the integrity of the kit before use.
- Use adequate protections (overall, gloves, glasses).
- Do not pipette by mouth.
- In case of contact with skin or eyes, thoroughly wash affected areas with plenty of water and seek medical advice.
- Material Safety Data Sheet is available upon request.
- Waste disposal: Respect legislation in force in the country.

All specimens should be handled as potentially infectious, in accordance with good laboratory practices using appropriate precautions. Respect legislation in force in the country.

REAGENT PREPARATION

Reagents are ready for use.



STABILITY AND STORAGE

Store at 18-25°C in well recapped original vial and away from light.

- **Standard (vial R2):** transfer requested quantity, well recap the vial and store at 18-25°C.
- Stored and used as described, reagents (vial R1 and R2) are stable, when free from contamination, until expiry date stated on the label of the kit.

Discard any reagent if cloudy or if absorbance at 340 nm > 0.500.

Don't use reagents after expiry date stated on the label of the kit.

SPECIMEN COLLECTION AND HANDLING (2)

Unhemolysed serum or plasma:

Separate from erythrocytes within 1 hour after collection.

- Phosphorus is stable in serum for:
 - ✓ several days at 2-8°C.
 - ✓ several months at -15°C.

24h Urines:

Collect in acid washed, detergent free container. Acidified urines (pH < 3 with concentrated hydrochloric acid) should be diluted (1+9) in demineralised water free from phosphorus before performing the assay.

- Phosphorus is stable in acidified urines for:
 - ✓ 6 months.

INTERFERENCES (3)

1. To avoid contamination with environmental phosphorus, it is recommended to use disposal glass or plastic ware, carefully cleaned material with hydrochloric acid 0.1 N and well rinsed with distilled water.

2. Gross hemolysis, lipemia or icterus may cause falsely elevated results unless a specimen blank is used.

For a more comprehensive review of factors affecting this assay refer to the publication of Young D.S.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Basic medical analysis laboratory equipment.
2. Normal and pathological control sera.

CALIBRATION (6)

- Standard (vial R2) provided in the kit or BIOLABO Multicalibrator REF 95015 traceable to SRM 3186a.
- Or any calibrator traceable to reference method or material.

The calibration frequency depends on proper instrument functions and on preservation of the reagent.

It is recommended to calibrate in the following cases :

1. When changing batch of reagent.
2. After maintenance operations on the instrument.
3. When control values are out of range, even after using a new vial of fresh serum.



Manufacturer



Use by



In vitro diagnostic



Temperature limitation



Catalogue number



See insert



Batch number



Store away from light



sufficient for



dilute with

QUALITY CONTROL

- BIOLABO EXATROL-N Level I **REF** 95010.
- BIOLABO EXATROL-P Level II **REF** 95011.
- Other assayed control sera referring to the same method.
- External quality control program.

It is recommended to control in the following cases:

- At least once a run.
- At least once within 24 hours.
- When changing vial of reagent.
- After maintenance operations on the instrument.

If control is out of range, apply following actions:

1. Repeat the test with the same control.
2. If control is still out of range, prepare a fresh control serum and repeat the test.
3. If control is still out of range, use a new vial of calibrator or a fresh calibrator and repeat the test.
4. If control is still out of range, calibrate with a new vial of reagent.
5. If control is still out of range, please contact BIOLABO technical support or your local Agent.

EXPECTED VALUES (2)

In serum or plasma

Phosphorus	(mg/dL)	(mmol/L)
In cord	3.7-8.1	[1.20-2.62]
Premature	5.4-10.9	[1.74-3.52]
0-10 days	4.5-9.0	[1.45-2.91]
10 days-24 months	4.5-6.7	[1.45-2.16]
24 months-12 years	4.5-5.5	[1.45-1.78]
12-60 years	2.7-4.5	[0.87-1.45]
> 60 years, male	2.3-3.7	[0.74-1.20]
> 60 years, female	2.8-4.1	[0.90-1.32]

In 24 h Urines

Phosphorus	(g/24 h)	(mmol/24 h)
Constant daily diet (*)	< 1	[< 32.3]
Non-restricted diet	0.4-1.3	[12.9-42.0]

Each laboratory should establish its own normal ranges for the population that it serves.

(*)Constant daily diet: 0.9 at 1.5 g P (29-48 mmol P) and 10 mg Ca/Kg (0.25 mmol Ca/Kg).

PERFORMANCES CHARACTERISTICS

Within run n = 20	Medium level	High level
Mean mg/dL	3.83	6.27
S.D. mg/dL	0.07	0.04
C.V. %	1.8	0.6

Between run n = 20	Medium level	High level
Mean mg/dL	3.62	7.53
S.D. mg/dL	0.12	0.16
C.V. %	3.3	2.1

Detection limit: approximately 0.7 mg/dL.

Sensitivity for 5 mg/dL: approximately 0.500 Abs at 340 nm.

Comparison study with commercially available method:

$$y = 0.990x + 0.0504 \quad r = 0.9948$$

LINEARITY

The assay is linear up to 10 mg/dL (3.22 mmol/L).

Above, dilute specimen with demineralised water free from phosphorus and re-assay taking into account dilution factor to calculate the result. Linearity limit depends on specimen/reagent ratio.

MANUAL PROCEDURE

Pipette into well identified test tubes:				
	Blank	Specimen blank	Standard	Assay
Reagent	1 mL		1 mL	1 mL
Saline Solution		1 mL		
Demineralised water	20 µL			
Standard			20 µL	
Specimen		20 µL		20 µL

Mix well. Incubate for 2 minutes at room temperature.
Use a 1 cm path length cuvette and read Standard and assays absorbance at 340 nm (334-366) against Reagent blank.
Read Specimen blank against saline solution.

Notes:

- ✓ Specific procedures are available upon request for automated instruments. Please contact BIOLABO technical support.
- ✓ Specimen blanking: for cloudy or icteric serum and for bovine control sera and calibrator which are usually very coloured.
- ✓ Target values of control sera and calibrator may have been obtained with or without specimen blanking.

CALCULATION

Calculate the result as follows:

Serum or plasma:

$$\text{Result} = \frac{\text{Abs (assay)} - \text{Abs (Specimen blank)}}{\text{Abs (Standard)}} \times \text{Standard concentration}$$

Urine diluted 1+9: Multiply the result by 10 (dilution factor).

REFERENCES

- (1) TIETZ N.W. *Text book of clinical chemistry*, 3rd Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p. 1406-1457.
- (2) *Clinical Guide to Laboratory Test*, 4th Ed., N.W. TIETZ (2006) p.852-855
- (3) YOUNG D.S., *Effect of Drugs on Clinical laboratory Tests*, 4th Ed. (1995) p.3-456 to 3-462
- (4) DALY J. A. and ERTINGSHAUSEN G., *Clin. Chem., Direct method for inorganic phosphate determination*, (1972), 18, p.263-265
- (5) GAMST O.K., TRY K., *Scand. J. Clin. Lab. Invest.* (1980), 40, p.483-486
- (6) SRM: Standard Reference Material®