



BIOLABO
www.biolabo.fr

MANUFACTURER:
BIOLABO SAS,

Les Hautes Rives
02160, Maizy, France

CHOLESTEROL

CHOD PAP Method

Reagent for quantitative determination of Total Cholesterol in human serum or plasma

REF	80106	R1	2 x 100 mL	R2	2 x 100 mL	R3	1 x 5 mL
REF	87356	R1	10 x 100 mL	R2	10 x 100 mL	R3	1 x 5 mL
REF	87656	R1	6 x 500 mL	R2	6 x 500 mL	R3	1 x 10 mL

TECHNICAL SUPPORT AND ORDERS

Tel : (33) 03 23 25 15 50

Fax: (33) 03 23 256 256



IVD IN VITRO DIAGNOSTIC USE

CLINICAL SIGNIFICANCE (1) (2)


Total cholesterol assay, associated to assays of others lipids in serum is used in the diagnosis of hyperlipidemia. Increased levels are also seen in hepatic and thyroid disorders.

Total cholesterol assay associated to triglycerides, HDL-Cholesterol and LDL-Cholesterol determination is useful in the prediction of coronary heart diseases.

So this assay is used in the diagnosis and treatment of atherosclerotic diseases. Hypercholesterolemia can also be observed in certain cases of diabetes. Secondary disorders that elevate cholesterol levels, should be ruled prior to initiating therapy with cholesterol-lowering drugs.

PRINCIPLE (4)

Enzymatic method described by Allain and al., which reaction scheme is as follows:



REAGENTS COMPOSITION

Vial R1		BUFFER	
Phosphate buffer	100	mmol/L	
Chloro-4-phenol	5	mmol/L	
Sodium Cholate	2.3	mmol/L	
Triton x 100	1.5	mmol/L	
Preservative			
Vial R2		ENZYMES	
Cholesterol oxydase (CO)	≥ 100	IU/L	
Cholesterol esterase (CE)	≥ 170	IU/L	
Peroxydase (POD)	≥ 1200	IU/L	
4 - Amino – antipyrine (PAP)	0.25	mmol/L	
PEG 6000	167	µmol/L	

Vial R3		STANDARD	
Cholesterol	200	mg/dL	(5.17 mmol/L)

SAFETY CAUTIONS

BIOLABO reagents are designated for professional, in vitro diagnostic use.

- Verify the integrity of the contents before use.
 - Use adequate protections (overall, gloves, glasses).
 - Do not pipette by mouth.
 - In case of contact with skin and eyes, thoroughly wash affected areas with plenty of water and seek medical advise.
 - Reagents contain sodium azide (concentration < 0.1%) which may react with copper and lead plumbing. Flush with plenty of water when disposing.
 - Material Safety Data Sheet is available upon request.
 - Waste disposal: Respect legislation in force in the country.
- All specimens should be handled as potentially infectious, in accordance with good laboratory practices using appropriate precautions. Respect legislation in force in the country.

REAGENTS PREPARATION

Add promptly the content of vial R2 (Enzymes), into vial R1 (Buffer).

Mix gently until complete dissolution (approximately 2 minutes).

Vial R2: If appropriate, use a non-sharp instrument to remove aluminium cap.

STABILITY AND STORAGE

Store at 2-8°C, well recap in the original vial and away from light.

- Standard (vial R3): Transfer the requested quantity, recap and store at 2-8°C.
- Reagents are stable until expiry date stated on the label of the kit when stored and used as described in the insert.
- Working reagent is stable at least for 2 years.
- Discard any reagent if cloudy or if reagent blank at 500 nm > 0.400.
- Don't use working reagent after expiry date stated on the label of the Kit.

SPECIMEN COLLECTION AND HANDLING (2)

Serum or plasma (Heparin or EDTA).

Do not use oxalate, fluoride or citrate. Collect on fasting patient. Separate serum from cells within 2 hours.

Cholesterol is stable in the specimen for:

- 5-7 days at 2-8°C
- 3 months at -20°C
- Many years at -70°C.
- Avoid repeated freezing and thawing

INTERFERENCES (2) (3) (5)

- Ascorbic acid:** Negative interference above 5 mg/dL.
Haemoglobin: Positive interference above 33 mg/dL.
Bilirubin: Negative interference above 8.6 mg/dL.
Serum blank may increase this interference.
Lipemia: Low interference limited by CE lipase activity.

Enzymatic methods increase analytic specificity. CO also reacts with 3 β -hydroxycholesterols (insignificant quantity in human serum – i.e. DHEA, pregnenolone).

For a more comprehensive review of factors affecting this assay refer to the publication of Young D.S or N. W. Tietz.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Basic medical analysis laboratory equipment.
2. Normal and pathological control sera.

CALIBRATION (6)

- Kit standard (vial R3) or BIOLABO Multicalibrator [REF] 95015
- Traceable to SRM 909b.
- Or any calibrator traceable to reference method or material.

The calibration frequency depends on proper instrument functions and on the preservation of the reagent.

It is recommended to calibrate in the following cases:

1. When changing vial of reagent.
2. After maintenance operations on the instrument.
3. When control values are out of range, even after using a new vial of fresh serum.

QUALITY CONTROL

- BIOLABO EXATROL-N Level I [REF] 95010.
- BIOLABO EXATROL-P Level II [REF] 95011.
- Other assayed control sera referring to the same method.
- External quality control program.

It is recommended to control in the following cases:

- At least once a run.
- At least once within 24 hours.
- When changing vial of reagent.
- After maintenance operations on the instrument.

If control is out of range, apply following actions:

1. Repeat the test with the same control.
2. If control is still out of range, prepare a fresh control serum and repeat the test.
3. If control is still out of range, use a new vial of calibrator or a fresh calibrator and repeat the test.
4. If control is still out of range, calibrate with a new vial of reagent.
5. If control is still out of range, please contact BIOLABO technical support or your local Agent.

EXPECTED VALUES (2)

Values for adults, estimated in term of risk for atherosclerotic diseases:

Total cholesterol	mg/dL	[mmol/L]
Recommended values	< 200	[< 5.18]
Low risk	200-239	[5.18-6.19]
High risk	\geq 240	[\geq 6.22]

Each laboratory should establish its own normal ranges for the population that it serves.

PERFORMANCE CHARACTERISTICS (4)

Within run N = 30	Normal level	High level	Between run: N = 33	Normal level	High level
Mean mg/dL	149	217	Mean mg/dL	125	261
S.D. mg/dL	1.25	1.89	S.D. mg/dL	1.04	2.06
C.V. %	0.84	0.87	C.V. %	0.83	0.79

Detection limit: approximately 1 mg/dL

Sensitivity for 100 mg/dL: 0.235 ± 0.035

Comparison study with commercially available reagent:

$$y = 0.957x + 6.4 \quad r = 0.9904$$

LINEARITY

The reaction is linear up to at least 500 mg/dL (12.9 mmol/L). Above, dilute the specimen with saline solution and re-assay taking into account dilution factor. Linearity limit depends on specimen/reagent ratio.

MANUAL PROCEDURE

Let stand reagent and specimens at room temperature.

Pipette into well identified test tubes:	Blank	Standard	Assay
Reagent	1 mL	1 mL	1 mL
Demineralised water	10 μ L		
Standard		10 μ L	
Specimen			10 μ L

Mix. Let stand for 5 minutes at 37°C or 10 minutes at room temperature. Record absorbances at 500 nm (480-520) against reagent blank. Colour is stable for 1 hour.

Note: Specific procedures are available upon request for automated instruments. Please contact BIOLABO technical support.

CALCULATION

Calculate the result as follows:

$$\text{Result} = \frac{\text{Abs (Assay)}}{\text{Abs (Standard)}} \times \text{Standard concentration}$$

REFERENCES

- (1) TIEZ N.W. *Text book of clinical chemistry*, 3rd Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p. 826-835.
- (2) *Clinical Guide to Laboratory Test*, 4th Ed., N.W. Tietz (2006) p. 244-249.
- (3) YOUNG D.S., *Effect of Drugs on Clinical laboratory Tests*, 4th Ed. (1995) p. 3-143 to 3-164
- (4) Allain C. C. et al., *Clin. Chem.* (1974), 20/4, p.470-475
- (5) Allan C., Deacon et Peter J. G. Dawson, *Clin. Chem.* (1979) 25/6, p.976-984
- (6) SRM: Standard Reference Material®