



BIOLABO
www.biolabo.fr

MANUFACTURER:
BIOLABO SAS,

Les Hautes Rives
02160, Maizy, France

TOTAL AND DIRECT BILIRUBIN

Sulfanilic Acid Method

Reagent for quantitative determination of total bilirubin (DMSO as accelerator)
or direct bilirubin in human serum and plasma

| | | | |
|-----|--------|---------------|------------------|
| REF | 80403: | R1 1 x 200 mL | Total bilirubin |
| | | R2 1 x 200 mL | Direct bilirubin |
| | | R3 1 x 40 mL | Nitrite Solution |
| REF | 80443: | R1 2 x 200 mL | Total bilirubin |
| | | R3 1 x 40 mL | Nitrite Solution |
| REF | 80553: | R2 2 x 200 mL | Direct bilirubin |
| | | R3 1 x 40 mL | Nitrite Solution |

TECHNICAL SUPPORT AND ORDERS

Tel: (33) 03 23 25 15 50

Fax: (33) 03 23 256 256



IVD IN VITRO DIAGNOSTIC USE

CLINICAL SIGNIFICANCE (1) (6)

At least four distinct bilirubin species exist in the serum: Direct-reacting bilirubin (**DB**) consists of mono and diconjugated bilirubin (β and γ -Bilirubin) and the δ -fraction which is bilirubin tightly bound to albumin; unconjugated α -bilirubin which is water soluble and bound to albumin. **Total bilirubin (TB)** is the sum of these different species.

There are icteruses in which unconjugated bilirubin predominates (hemolytic icteruses, Biermer disease, Thalassemia...); icteruses in which conjugated bilirubin predominates (extra or intra-hepatic bile ducts obstruction, viral hepatitis...); finally, icteruses in which both species of bilirubin are present without any predominance (cirrhosis, Dubin-Johnson disease).

PRINCIPLE (4) (5)

Reaction between bilirubin and diazotised sulfanilic acid which leads to a compound, the azobilirubin, coloured in very acid or basic medium.

Malloy-Evelyn principle modified by Walters and al: in an aqueous solution, only DB reacts. To enable the assay of TB, it is necessary to break the link between unconjugated bilirubin and albumin. This step is done by adding dimethylsulfoxide (DMSO).

The absorbance of azobilirubin thus produced is proportional to the concentration of bilirubin and is measured at 550 nm (530-580).

REAGENTS COMPOSITION

| Vial R1 | TOTAL BILIRUBIN | |
|-------------------|------------------|--------|
| Sulfanilic acid | 30 | mmol/L |
| DMSO | 7 | mol/L |
| Hydrochloric acid | 130 | mmol/L |
| Vial R2 | DIRECT BILIRUBIN | |
| Sulfanilic acid | 30 | mmol/L |
| Hydrochloric acid | 130 | mmol/L |
| Vial R3 | NITRITE SOLUTION | |
| Sodium Nitrite | 0.74 | mmol/L |

SAFETY CAUTIONS

BIOLABO reagents are designated for professional, in vitro diagnostic use.

- Verify the integrity of the contents before use.
- Use adequate protections (overall, gloves, glasses).
- Do not pipette by mouth.
- In case of contact with skin and eyes, thoroughly wash affected areas with plenty of water and seek medical advice.
- Material Safety Data Sheet is available upon request.
- Waste disposal: Respect legislation in force in the country.

All specimens should be handled as potentially infectious, in accordance with good laboratory practices using appropriate precautions. Respect legislation in force in the country.



REAGENTS PREPARATION

Reagents are ready for use.
In case of important runs or automated system one can prepare single vial working reagent as follows: R1 or R2 (20 volumes)+R3 (1 volume).

STABILITY AND STORAGE

Store at 2-8°C, in well capped original vial and away from light.

- When used and stored as described in the insert, reagents (vial R1, R2, R3) are stable, without contamination, until expiry date stated on the label.
- Working reagent TB is stable for 2 days at 2-8°C.
- Working reagent DB is stable for 7 days at 2-8°C.

Discard any reagent if cloudy or if absorbance at 550 nm > 0.100.

Don't use working reagents after expiry date stated on the label of the kit.

SPECIMEN COLLECTION AND HANDLING (2) (7)

Unhemolysed serum or plasma.

Bilirubin is photolabile. Store the specimen away from light.

- Stability in the specimen: 4 to 7 days at 2-8°C.
2 days at room temperature.

Icteric or pediatric specimens: see § **MANUAL PROCEDURE**.

INTERFERENCES (3)

Hemoglobin: under-evaluation above 100 μ mol/L (160 mg/dL) of haemoglobin.

Turbidity: No significant interference with TB.
No significant interference with DB up to triglycerides concentration equivalent to 4.6 mmol/L.

The bilirubin diazo reaction is temperature sensitive and should be carried out at a constant temperature.

For a more comprehensive review of factors affecting this assay refer to the publication of Young D.S.

MATERIAL REQUIRED BUT NOT PROVIDED

1. Basic medical analysis laboratory equipment.
2. Normal and pathological control sera.

CALIBRATION (8)

1. Use the experimental factor (see § **CALCULATION**)
2. Or a calibrator traceable to a reference method or material, BIOLABO Multicalibrator REF 95015 traceable to SRM 916a.

The calibration frequency depends on proper instrument functions and on preservation of reagents.

It is recommended to calibrate in the following cases:

1. When using a new batch of reagent.
2. After maintenance operations on the instrument.
3. When control values obtained are out of ranges, even after using a new vial of fresh serum



Manufacturer



Use by



In vitro diagnostic



Temperature limitation



Catalogue number



See insert



Batch number



Store away from light



sufficient for



dilute with

QUALITY CONTROL

- BIOLABO EXATROL-N Level I [REF] 95010.
- BIOLABO EXATROL-P Level II [REF] 95011.
- BIOLABO PAEDIATRIC CONTROL [REF] 95403
- Other assayed control sera referring to the same method and selected procedure.
- External quality control program.

It is recommended to control in the following cases:

- At least once a run.
- At least once within 24 hours.
- When changing vial of reagent.
- After maintenance operations on the instrument.

If control is out of range, apply following actions:

1. Repeat the test with the same control.
 2. If control is still out of range, prepare a fresh control serum and repeat the test.
 3. **With factor:** Verify analysis parameters (Wavelength, temperature, specimen/reagent ratio, time counting, calibration factor).
 4. Use a new vial of reagent and repeat the test.
 5. **With a calibrator:** If control is still out of range, use a new vial of calibrator or a fresh calibrator and repeat the test.
 6. If control is still out of range, calibrate again with a new vial of reagent and repeat the test.
1. If control is still out of range, please contact BIOLABO technical support or your local Agent.

EXPECTED VALUES (2)

| Total Bilirubin | mg/dL | | [µmol/L] | |
|-----------------|-----------|-----------|-----------|-----------|
| Newborn | Premature | Full-term | Premature | Full-term |
| In cord | < 2.0 | < 2.0 | [< 34] | [< 34] |
| 0-1 day | < 8.0 | 1.4-8.7 | [< 137] | [24-149] |
| 1-2 days | < 12.0 | 3.4-11.5 | [< 205] | [58-197] |
| 3-5 days | < 16.0 | 1.5-12.0 | [< 274] | [26-205] |

| Adult (and child > 5 days) | Total Bilirubin | | Direct Bilirubin | |
|----------------------------------|-----------------|----------|------------------|----------|
| | mg/dL | [µmol/L] | mg/dl | [µmol/L] |
| >5 days-60 years | 0.3-1.2 | [5-21] | < 0.2 | [< 3.4] |
| 60-90 years | 0.2-1.1 | [3-19] | < 0.2 | [< 3.4] |
| > 90 years | 0.2-0.9 | [3-15] | < 0.2 | [< 3.4] |

Each laboratory should establish its own normal ranges for the population that it serves.

LINEARITY

Procedure n°1: up to 20 mg/dL (342 µmol/L).

Above, **do not dilute specimen:** perform procedure n°2.

Procedure n°2: up to 100 mg/dL (1710 µmol/L)

Pediatric specimen: perform procedure n°2

PERFORMANCES CHARACTERISTICS (PROCEDURE N°1)

TOTAL BILIRUBIN

| Within run: N = 23 | Normal level | High level | Between run N = 20 | Normal level | High level |
|-----------------------|--------------|------------|-----------------------|--------------|------------|
| | Mean mg/dL | 0.68 | | 4.13 | Mean mg/dL |
| S.D. mg/dL: | 0.02 | 0.104 | S.D. mg/dL: | 0.022 | 0.07 |
| C.V. %: | 2.94 | 2.52 | C.V. %: | 3.27 | 1.78 |

DIRECT BILIRUBIN

| Within run: N = 20 | Medium level | High level | Between run N = 20 | Medium level | High level |
|-----------------------|--------------|------------|-----------------------|--------------|------------|
| | Mean mg/dL | 1.15 | | 2.79 | Mean mg/dL |
| S.D. mg/dL: | 0.022 | 0.015 | S.D. mg/dL: | 0.028 | 0.09 |
| C.V. %: | 1.94 | 0.53 | C.V. %: | 2.6 | 3.3 |

Detection limit: TB: approximately 0.13 mg/dL

DB: approximately 0.18 mg/dL

Sensitivity for 1 mg/dL: 88 mAbs at 550 nm.

Comparison study with commercially available reagent:

TB: $y = 1.0145x + 0.00513$ $r = 0.9976$

DB: $y = 1.0002x - 0.00796$ $r = 0.9972$

MANUAL PROCEDURE

Let stand reagents and specimens at room temperature.

Procedure n°1:

| Pipette into well identified test tubes: | TOTAL BILIRUBIN | | DIRECT BILIRUBIN | |
|--|-----------------|--------|------------------|--------|
| | Blank | Assay | Blank | Assay |
| Reagent R1 | 1 mL | 1 mL | | |
| Reagent R2 | | | 1 mL | 1 mL |
| Distilled water | 50 µL | | 50 µL | |
| Reagent R3 (Nitrite) | | 50 µL | | 50 µL |
| Mix | | | | |
| Specimen | 100 µL | 100 µL | 100 µL | 100 µL |

Mix well and start a timer when adding specimen.
Read absorbances at 550 nm (530-580) against blanks
TB: reading after ≥ 3 minutes at 37°C or ≥ 5 minutes at room temperature.
DB: reading at **exactly** 3 minutes at 37°C or 5 minutes at room temperature.

Procedure n°2: Icteric or Pediatric Specimens

| Pipette into well identified test tubes | TOTAL BILIRUBIN | | DIRECT BILIRUBIN | |
|---|-----------------|-------|------------------|-------|
| | Blank | Assay | Blank | Assay |
| Reagent R1 | 1 mL | 1 mL | | |
| Reagent R2 | | | 1 mL | 1 mL |
| Distilled water | 50 µL | | 50 µL | |
| Reagent R3 (Nitrite) | | 50 µL | | 50 µL |
| Mix | | | | |
| Specimen | 20 µL | 20 µL | 20 µL | 20 µL |

Well mix and start a timer when adding specimen.
Read absorbance at 550 nm (530-580) against blanks
TB: reading after ≥ 3 minutes at 37°C or ≥ 5 minutes at room temperature.
DB: reading at **exactly** 3 minutes at 37°C or 5 minutes at room temperature.

Notes:

1. Zero on distilled water and **drain well** the cuvette. Read first all blanks of one run then all the assays, well draining cuvette between each tubes. But **DO NOT RINSE WITH WATER** as it could produce streaks on the cuvette and lead to false results.
2. Specific procedures are available upon request for automated instruments. Please contact BIOLABO technical support.

CALCULATION

Calculate the result as follows:

With calibrator (Procedure n°1 only):

$$\text{Result} = \frac{\text{Abs (Assay - Blank) Specimen}}{\text{Abs (Assay - Blank) Calibrator}} \times \text{calibrator concentration}$$

With factor:

Procedure n°1: $\text{mg/dL} = [\text{Abs. assay} - \text{Abs. Blank}] \times 11.4^*$

$\mu\text{mol/L} = [\text{Abs. assay} - \text{Abs. Blank}] \times 195^*$

Procedure n°2: $\text{mg/dL} = [\text{Abs. assay} - \text{Abs. Blank}] \times 53.0^*$
 $\mu\text{mol/L} = [\text{Abs. assay} - \text{Abs. Blank}] \times 906^*$

*This factors should be used as a guide only and may vary with instrument and the batch of reagent used. It is recommended to verify with elevated control serum

(Procedure n°1: Use BIOLABO EXATROL-P, Procedure n°2: Use BIOLABO PAEDIATRIC CONTROL).

REFERENCES

- (1) TIETZ N.W. *Text book of clinical chemistry*, 3rd Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p. 1133-1137.
- (2) *Clinical Guide to Laboratory Test*, 4th Ed., N.W. TIETZ (2006) p. 172-177.
- (3) YOUNG D.S., *Effect of Drugs on Clinical laboratory Tests*, 4th Ed. (1995) p. 3-90 to 3-110
- (4) MALLOY H.T., EVELYN K., *J Biol. Chem.*(1937), 119, p.481-490
- (5) WALTERS M, GERARDE H, *Microchem J* (1970) 15 p.231-243
- (6) BERNARD S., *Biochimie clinique*, 2^{ème} éd. Maloine, (1989), p.127-129 et p.280-282
- (7) Henry RJ, *Clin Chem: Principles and technics*. Harper and Row. p.592(1965).
- (8) SRM: Standard Reference Material®