



**BIOLABO**  
www.biolabo.fr

**MANUFACTURER:**  
**BIOLABO SAS,**

Les Hautes Rives  
02160, Maizy, France

# Lyophilised G6-PDH U.V. Kinetic Method

Reagent for quantitative determination of Glucose-6-phosphate dehydrogenase activity (G6-PDH)  
[ EC 1.1.1.49 ] in human serum, plasma or erythrocytes

REF 97099 R1 20 x 3 mL R2 3 x 1 mL R3 1 x 70 mL R4 1 x 20 mL

## TECHNICAL SUPPORT AND ORDERS

Tel: (33) 03 23 25 15 50

Fax: (33) 03 23 256 256



IVD IN VITRO DIAGNOSTIC USE

## CLINICAL SIGNIFICANCE (1) (2)


Glucose-6-phosphate dehydrogenase (G6-PDH) is an enzyme of glycolysis pathway in erythrocytes.

Different variants of G6-PDH defined by electrophoretic and kinetic criteria are of high frequency and have clinical consequences in different ethnic groups.

A G6-PDH deficiency may cause haemolytic anemia following ingestion of 8-amino-quinoline antimalarials, nalidixic acid, nitrofurantoin, phenacetin, large doses of vitamin C and some sulfonamides and sulfones. Deficiency may be the cause of haemolytic disease of newborns in Asians and Mediterraneans.

## PRINCIPLE (4) (5)

Reaction scheme (method of Beutler and al.) is as follows:



The rate of increase in NADPH concentration measured at 340 nm is proportional to the G-6-PDH activity in the specimen.

## REAGENTS

<b>Vial R1</b>	<b>COENZYME</b> Freeze-dried
NADP <sup>+</sup>	310 mmo/L
<b>Vial R2</b>	<b>SUBSTRATE</b> Freeze-dried
Glucose-6-phosphate	0.6 mmol/L
<b>Vial R3</b>	<b>BUFFER</b>
Tris Buffer pH 8.0	100 mmol/L
MgCl <sub>2</sub>	10 mmol/L
EDTA	0.5 mmol/L
Preservative	
<b>flacon R4</b>	<b>HEMOLYSING SOLUTION</b>
Digitonin	0.2 g/L

## SAFETY CAUTIONS

BIOLABO reagents are designated for professional, in vitro diagnostic use.

- Use adequate protections (overall, gloves, glasses).
- Do not pipette by mouth.
- In case of contact with skin or eyes, thoroughly wash affected areas with plenty of water and seek medical advice.
- Reagents contain sodium azide (concentration < 0.1%) which may react with copper and lead plumbing. Flush with plenty of water when disposing.
- Material Safety Data Sheet is available upon request.
- Waste disposal: Respect legislation in force in the country.

All specimens should be handled as potentially infectious, in accordance with good laboratory practices using appropriate precautions. Respect legislation in force in the country.

## REAGENTS PREPARATION

Add promptly 3 mL of buffer (vial R3) to the contents of vial R1 (Coenzyme) and 1 mL of buffer (vial R3) to the contents of vial R2 (Substrate).

Mix gently and wait for complete dissolution before using reagents (approximately 2 minutes).

Hemolysing solution (vial R4) is ready for use.

## STABILITY AND STORAGE

Store at 2-8°C, well recapped in the original vial and away from light.

- Stored and used as described in the insert, unopened reagents are stable until expiry date stated on the label.
- After reconstitution working reagents (vial R1 and R2) are stable at least for 6 months when free from contamination.
- Don't use working reagents after expiry date stated on the label. Discard any reagent if cloudy.

## SPECIMEN COLLECTION AND HANDLING

Unhemolysed serum or plasma (heparinised or EDTA). Avoid oxalate and fluoride.

Erythrocytes: Use an hemolysate of whole blood prepared as follows: (Detailed Procedure available upon request)

1. Determine the concentration of haemoglobin (Hb in g/dL).
2. Wash 3 times 0,2 mL of homogenised blood with 2 mL of saline solution 0.9 g/dL. Centrifuge between each washing and eliminate the supernatant (avoiding elimination of erythrocytes).
3. After the last washing, suspend the washed erythrocytes in 0,9 mL of haemolysing solution (vial R4)
4. Let stand for 15 minutes at 2-8°C and centrifuge again. Use the supernatant (hemolysate) within 1 hour.

## INTERFERENCES (2) (3)

This assay reflects also the 6-Phosphogluconate dehydrogenase (6-PGD) activity which generates one molecule of NADPH for one molecule of 6-Phosphogluconate formed.

G6PDH (UI/L) at 37°C in the specimen	Interfering	Results
169 UI/L	Ascorbic acid	Do not interfere up to 25 mg/dL
1756,7 UI/L	Turbidity	Do not interfere up to 0,2 %
1846,7 UI/L	Glucose	Do not interfere up to 1100 mg/dL
2003,3 UI/L	Bilirubin	Negative interference above 7,5 mg/dL

For a more comprehensive review of factors affecting this assay refer to the publication of Young D.S.

## MATERIAL REQUIRED BUT NOT PROVIDED

1. Basic medical analysis laboratory equipment.
2. Normal and pathological control sera or hemolysate



Manufacturer



Use by



In vitro diagnostic




Temperature limitation



Catalogue number



See insert



Batch number



Store away from light



sufficient for



dilute with

## CALIBRATION

Results will depend on the accuracy of the instrument calibration, the time counting, the respect of reagent/specimen ratio and the temperature control.

It is recommended to use the theoretical calibration factor (§ **CALCULATION**) or a seric calibrator traceable to a reference method or a reference hemolysate.

## QUALITY CONTROL

- **REF** 95089 Normal Control G6-PDH
  - **REF** 95289 Deficient Control G6-PDH
  - Any assayed control sera or a whole blood specimen referring to the same method.
  - External quality control program.
- It is recommended to control in the following cases:
- At least once a run.
  - At least once within 24 hours.
  - When changing vial of reagent.
  - After maintenance operations on the instrument.

If control is out of range, apply following actions:

1. Repeat the test with the same control.
2. If control is still out of range, prepare a fresh control and repeat the test.
3. If control is still out of range, verify analysis parameters: Wavelength, temperature, specimen/reagent ratio, time counting, calibration factor.
4. If control is still out of range, use a new vial of reagent and reassay
5. If control is still out of range, please contact BIOLABO technical support or your local Agent.

## EXPECTED VALUES (2)

### In Erythrocytes at 37°C

Conventional Units	International Units
IU/g of Hb: 12.1 ± 2.09	MIU/molHb: [0.78 ± 0.13]
IU/10 <sup>12</sup> Erythrocytes: 351 ± 60.6	nIU/Erythrocytes: [0.35 ± 0.06]
IU/mL of Erythrocytes: 4.11 ± 0.71	KIU/L of Erythrocytes: [4.11 ± 0.71]

Each laboratory should establish its own normal ranges for the population that it serves..

### In Serum at 37°C

Normal value: No detectable G6-PDH activity.

## PERFORMANCE CHARACTERISTICS

Studies performed using hemolyzate as specimen (Cobas Mira).

Within run n=20	Low level	Normal level	High level
Mean (IU/L)	478	925	2280
S.D. (IU/L)	13.3	33.2	31.2
C.V. %	2.8	3.6	1.4
Critères CV%	< 4.5%	< 4.5%	< 3.8%

Between run	Low level (n=15)	Normal level (n=15)	High level (n=30)
Mean (IU/L)	276	535	918
C.V. %	5.8	4.8	3.9
S.D (IU/L)	16.0	25.5	35.5
C.V. %	5.8	4.8	3.9
Critères CV%	< 6%	< 6%	< 5%

Detection limits: approximately 21 IU/L.

Sensitivity: approx. 0.020 ΔAbs/min /1000 IU.L<sup>-1</sup> of blood

Comparison study with commercially available reagent (UV Kinetic method):

102 hemolysates within 110 and 1500 IU/L have been evaluated with both reagents on Cobas Mira:

$$y = 0.9992x - 6.5933$$

$$r = 0.9874$$

## LINEARITY

Above 4000 UI/L (0.080 ΔAbs/mn), dilute specimen (serum, plasma, hemolysate) with saline solution 0.9 g/dL and reassay taking into account dilution factor to calculate the result. Linearity limits depend on the specimen/reagent ratio.

## MANUAL PROCEDURE

Let stand reagents and specimens at room temperature.

Pipette into thermostated cuvette at 37°C:	Serum Assay	Hemolysate Assay(1)
<b>Reagent R1</b>	2 mL	3 mL
<b>Serum</b>	1 mL	
<b>Hemolysate</b>		50 µl

Mix and incubate for 5 minutes at 37°C (30°C, 25°C)

Reagent R2	100 µl	100 µl
------------	--------	--------

Mix and record initial absorbance at 340 nm after 30 seconds. Record the absorbance again every minutes during 3 minutes.  
Calculate absorbance change per minute (ΔAbs/min).  
If activity is lower, timing of measurement may be prolonged.

Notes (1):

1-Preparation of the hemolysate: see § **SPECIMEN COLLECTION AND HANDLING**.

2- Hemolysate assay may be realised directly in the vial R1

## CALCULATION (2)

Calculate G6-PDH activity as follows:

**Serum:** IU/L = (ΔAbs/min) x 492

**Erythrocytes:** IU/L of blood = (ΔAbs/min) x 50 000

**Result in units per gramme of Hemoglobin**

$$\text{IU/g Hb} = \frac{(\Delta\text{Abs/min} \times 5\,000)}{\text{Hb expressed in g/dL}}$$

Example: if ΔAbs/min = 0.030 and Hb = 14.5 g/dL

$$\text{IU/g Hb} = \frac{0.030 \times 5000}{14.5} = 10.3$$

**Result in units per 10<sup>12</sup> Erythrocytes (# Result in IU/g Hb x 29)**

$$\text{IU/10}^{12} \text{ Erythrocytes} = \frac{\Delta\text{Abs/min} \times 50\,000}{\text{Number of erythrocytes in } 10^{12}/\text{L}}$$

Example: if ΔAbs/min = 0.030 and erythrocytes number = 4.2. 10<sup>12</sup>/L

$$\text{IU/10}^{12} \text{ erythrocytes} = \frac{0.030 \times 50\,000}{4.2} = 357$$

**Result in units per mL of Erythrocytes (# Result in IU/g Hb x 0.34)**

$$\text{IU/ml of Erythrocytes} = \frac{(\Delta\text{Abs/min} \times 5\,000)}{\text{Hematocrit (\%)}}$$

## REFERENCES

- (1) TIETZ N.W. *Textbook of clinical chemistry*, 3<sup>rd</sup> Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p. 1645-1650.
- (2) *Clinical Guide to Laboratory Test*, 4<sup>th</sup> Ed., N.W. TIETZ (2006) p. 457-458.
- (3) YOUNG D.S., *Effect of Drugs on Clinical Laboratory Tests*, 4<sup>th</sup> Ed. (1995) p. 3-294
- (4) BEUTLER et al., *International committee for standardisation in Haematology: « Recommended Methods for Red Cell Enzyme Analysis » British Journal of Haematology*, (1977), 35, p.331-340.
- (5) BEUTLER E., *Red cell metabolism. A manual of biochemical methods* (3<sup>rd</sup> Ed.) Orlando, Grune et Stratton (1984), p.68-70