



BIOLABO
www.biolabo.fr
MANUFACTURER:
BIOLABO SAS,
Les Hautes Rives
02160, Maizy, France

CONTROL PLASMA Level 2

For Internal Quality Control in Haemostasis

REF 13962 R1 6 x 1 mL



Made In France

TECHNICAL SUPPORT AND ORDERS
Tel: (33) 03 23 25 15 50
Fax : (33) 03 23 256 256
support@biolabo.fr
Latest revision: www.biolabo.fr

I: corresponds to significant modifications

VALUES AND RANGES ⁽³⁾

CONTROL PLASMA Level 2

LOT _____

Make sure that the batch number stated on the label of the vial corresponds to the batch number indicated here above.
Fibrinogen is determined against secondary standard of SSC/ISTH Secondary Coagulation Standard NIBSC Code: SSCLOT4

		Semi-automated and full-automated methods, including: BIOLABO SOLEA 100, BEHNK THROMBOLYZER Series	
		Target value	Confidence Interval
BIO -TP LI:	INR		
Low ISI	Prothrombin level or PT (%)		
BIO-TP:	INR		
High ISI	Prothrombin level or PT (%)		
BIO-CK:	Activated Partial Thromboplastin Time (sec)		
BIO-SIL:	Activated Partial Thromboplastin Time (sec)		
BIO-FIBRI:	Fibrinogen (mg/dL)		

INTENDED USE

| This control plasma is designated for professional use in laboratory (manual or automated method).
It is used to monitor the reproducibility and accuracy of following analysis performed with BIOLABO Reagents:

- REF 13702, 13704 and 13712 : BIO-TP LI
- REF 13885, 13880 and 13881 : BIO-TP
- REF 13560 and 13570 : BIO-CK
- REF 13660 and 13670 : BIO-SIL
- REF 13450 and 13451: BIO-FIBRI

REAGENTS

REF 13962: **CONTROL PLASMA** Level 2



Freeze-dried, human citrated plasmas

Human Origin

SAFETY CAUTIONS ^{(1) (2)}

- Material Safety Data Sheet is available upon request.
- Each donor unit used to manufacture this product was tested and found non-reactive for HbsAg, antibody to Hepatitis C and antibody to HIV-1/HIV-2.
- However, no test method can offer complete assurance that infectious agents are absent. All specimens or reagents from biological origin should be handled as potentially infectious, in accordance with good laboratory practices using appropriate precautions.
- Waste disposal: Respect legislation in force in the country.

| Any serious incident that has occurred in connection with the device is notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is based.

MATERIAL REQUIRED BUT NOT PROVIDED

1. Basic medical analysis laboratory equipment
2. Control Plasmas REF 13961, REF 13963
3. Reference plasma REF 13970 or TP CAL SET REF 13965

REAGENTS PREPARATION

- Open the vial carefully and add exactly the volume of demineralised water **stated on the label**
- Wait for 15 minutes at room temperature.
- Gently agitate before use (avoid the formation of foam).

WARNING: DO NOT SHAKE. STORE AWAY FROM LIGHT

STABILITY AND STORAGE

Stored away from light, well capped in the original vial at 2-8°C, when stored and used as described, reagents are stable:

Unopened:

- Until expiry date stated on the label of the kit.

Once opened:

- R1 must be reconstituted immediately,
- | Once reconstituted, plasma is stable:
 - 10 hours at 2-25°C
 - 5 days at -20°C

Don't use reconstituted plasma after expiry date stated on the label

PROCEDURE

Run in accordance with the IFU of the reagent used § **CONTROL**

LIMITES

Factors which may interfere with the result are:

- Bacterial contamination.
- The volume measured to reconstitute the plasma.
- The setting of the instrument.
- Temperatures

REFERENCES

- (1) *Occupational Safety and Health Standards: Bloodborne pathogens (29CFR1910.1030) Federal Register July 1, (1998) ; 6, p.267-280*
- (2) *Directive du conseil de l'Europe (90/679/CEE) J. O. de la communauté européenne n°L374 du 31.12.1990, p.1-12*
- (3) *Section 5.6 of ISO 17511- Measurements of quantities in biological samples- metrological traceability of values assigned to calibrators and contro*

Manufacturer REF	Expiry date See Insert	In vitro diagnostic LOT	Storage temperature Store away from light	Dematerialized water Sufficient for	Biological risk → Dilute with
---------------------	---------------------------	----------------------------	--	--	----------------------------------