



BIOLABO
 www.biolabo.fr
MANUFACTURER:
BIOLABO SAS,
 Les Hautes Rives
 02160, Maizy, France

ASLO Standard Set

For calibration of quantitative determination
 of Antistreptolysin O (ASLO) by Turbidimetric Immunoassay

REF ASLO CALSET41
 R1 1 x 1 mL, R2 1 x 1 mL, R3 1 x 1 mL, R4 1 x 1 mL



Made in France

I: corresponds to significant modifications

TECHNICAL SUPPORT AND ORDERS

Tel : (33) 03 23 25 15 50

Fax: (33) 03 23 256 256

support@biolabo.fr

Latest revision : www.biolabo.fr

INTENDED USE

Standard Set for the preparation of reference curve for quantitative immunochemical determination of Antistreptolysin O in human serum. Suitable for manual procedure or automated instruments with BIOLABO reagents REF ASLO050E, REF ASLO620E.

REAGENTS

R1	ASLO CAL1	 Human origin
R2	ASLO CAL2	
R3	ASLO CAL3	
R4	ASLO CAL4	

4 vials of ASLO Standards (4 different levels)

Liquid plasmas supplemented with ASLO, diluted in Saline and stabilized.

SAFETY CAUTIONS (1) (2)

- BIOLABO reagents are designated for professional use in laboratory.
- Material Safety Data Sheet is available upon request.
 - Each donor unit used to manufacture this product was tested and found non-reactive for HbsAg, antibody to Hepatitis C and antibody to HIV-1/HIV-2.
 - However, no test method can offer complete assurance that infectious agents are absent. All specimens or reagents from biological origin should be handled as potentially infectious, in accordance with good laboratory practices using appropriate precautions.
 - Waste disposal: Respect legislation in force in the country.
- I Any serious incident that has occurred in connection with the device is notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is based.

REAGENTS PREPARATION

Ready for use

MATERIAL REQUIRED BUT NOT PROVIDED

- BIOLABO Reagents (§ INTENDED USE)
- REF ASLO CONT1 or REF ASLO CONT5: ASLO Control

QUALITY CONTROL

Verify the integrity of each vials and batch-specific values before use. Run in accordance with the IFU of the reagent used.



STABILITY AND STORAGE

Stored away from light, well cap in the original vial at 2-8°C, standards are stable when stored and used as described:

- Unopened,
- Until the expiry date stated on the label of the Kit.
- Once opened:
- Transfer requested quantity, well recap vials and store at 2-8°C.
 - Well recapped in the original vial, at least for 6 weeks when free from contamination.

Do not freeze

PROCEDURE

Run in accordance with the IFU of the reagent used.

CALIBRATION VALUES (3)

- Values are traceable to a reference material (WHO Standardisation)
- Batch-specific** values are indicated on the label of each vial

LIMITATIONS

Factors which may influence results are bacterial contamination, accuracy of reconstitution volume, respect of automated instrument procedure, temperature...

REFERENCES

- Occupational Safety and Health Standards ; Bloodborne pathogens (29CFR1910.1030) Federal Register July 1, (1998) ; 6, p.267-280
- Directive du conseil de l'Europe (90/679/CEE) J. O. de la communauté européenne n°L374 du 31.12.1990,p.1-12
- TIETZ N.W. Textbook of clinical chemistry, 3rd Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p.520

 Manufacturer	 Expiry date	 In vitro diagnostic	 Storage temperature	 Dematerialized water	 Biological risk
 Product Reference	 See Insert	 Batch number	 Store away from light	 Sufficient for	 Dilute with