



BIOLABO
www.biolabo.fr
MANUFACTURER:
BIOLABO SAS,
 Les Hautes Rives
 02160, Maizy, France

CRP Standard Super High

For calibration of quantitative determination
 of C reactive protein (CRP) by Turbidimetric immunoassay

REF CRP CALSH1 **R1** 1 x 1 mL



Made in France

TECHNICAL SUPPORT AND ORDERS

Tel : (33) 03 23 25 15 50

Fax: (33) 03 23 256 256

support@biolabo.fr

Latest revision : www.biolabo.fr

I: corresponds to significant modifications

INTENDED USE

Standard to be diluted for the calibration of quantitative immunochemical determination of CRP in human serum. Suitable for manual procedure or automated instruments with BIOLABO reagents **REF** CRP050E, **REF** CRP620E.

REAGENTS

R1 CRP Standard Super High



Human Origin

Liquid stabilized plasma (defibrinated and delipidated) supplemented with purified CRP.

SAFETY CAUTIONS (1) (2)

- BIOLABO reagents are designated for professional use in laboratory.
- Material Safety Data Sheet is available upon request.
 - Each donor unit used to manufacture this product was tested and found non-reactive for HbsAg, antibody to Hepatitis C and antibody to HIV-1/HIV-2.
 - However, no test method can offer complete assurance that infectious agents are absent. All specimens or reagents from biological origin should be handled as potentially infectious, in accordance with good laboratory practices using appropriate precautions.
 - Waste disposal: Respect legislation in force in the country.

I Any serious incident that has occurred in connection with the device is notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is based.

REAGENTS PREPARATION

Ready for use

MATERIAL REQUIRED BUT NOT PROVIDED

1. BIOLABO Reagents (§ INTENDED USE)
2. **REF** CRP CONT L1, CRP CONT L5: CRP Low Control
3. **REF** CRP CONT H1, CRP CONT H5: CRP High Control

QUALITY CONTROL

Verify the integrity of the vial and batch-specific value before use. Run in accordance with the IFU of the reagent used.

STABILITY AND STORAGE

Stored away from light, well cap in the original vial at 2-8°C, the standard is stable when stored and used as described:

Unopened,

- Until the expiry date stated on the label of the Kit.

Once opened:

- Transfer requested quantity, well recap vials and store at 2-8°C.
- Well recapped in the original vial, at least for 6 weeks when free from contamination.

Do not freeze

PROCEDURE

Run in accordance with the IFU of the reagent used

CALIBRATION VALUES (3)

- Standard value is traceable to a reference material (RPPHS/CRM470) from the International Federation of Clinical Chemistry (IFCC)
- **Batch-specific** value is indicated on the label of the vial

LIMITATIONS

Factors which may influence results are bacterial contamination, accuracy of reconstitution volume, respect of automated instrument procedure, temperature...

REFERENCES

- (1) *Occupational Safety and Health Standards ; Bloodborne pathogens (29CFR1910.1030) Federal Register July 1, (1998) ; 6, p.267-280*
- (2) *Directive du conseil de l'Europe (90/679/CEE) J. O. de la communauté européenne n°L374 du 31.12.1990,p.1-12*
- (3) *TIETZ N.W. Textbook of clinical chemistry, 3rd Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p.520*

Manufacturer REF Product Reference	Expiry date See Insert	In vitro diagnostic LOT Batch number	Storage temperature Store away from light	Dematerialized water Sufficient for	Biological risk → Dilute with
-------------------------------------------------	---------------------------	---------------------------------------------------	----------------------------------------------	----------------------------------------	-------------------------------------