



**BIOLABO**  
www.biolabo.fr

**MANUFACTURER:**  
**BIOLABO SAS,**

Les Hautes Rives  
02160, Maizy, France

# CRP Turbidimetric Immunoassay

Reagent for quantitative determination  
of C-Reactive Protein (CRP) in human serum.

REF CRP050E	R1 1 x 50 mL	R2 1 x 5 mL
REF CRP620E	R1 6 x 20 mL	R2 1 x 10 mL



Made in France

## TECHNICAL SUPPORT AND ORDERS

Tel: (33) 03 23 25 15 50

Fax: (33) 03 23 256 256

support@biolabo.fr

Latest revision : www.biolabo.fr

I: corresponds to significant modifications

## INTENDED USE

This quantitative test is to detect and measure C - reactive protein in human serum to assess the inflammatory status of the body. This test may be used with manual procedure on spectrophotometer or with Biochemistry Clinical Analyzer.

## GENERALITIES (1) (4)

C - reactive protein is one of the strongest acute phase reactants and help in non-specific host defence against infectious agents.

Its concentration increased after myocardial infarction, stress, trauma, infection, inflammation, surgery or neoplastic proliferation.

## PRINCIPLE (2) (3)

Turbidimetric Immunoassay (TIA): Photometric measurement of turbidity, corresponding to antigen-antibody reaction, by the end-point method at 340 nm.

## REAGENTS

R1	CRP TIA	Buffer
Phosphate buffered saline		pH 7.43
Polyethylene glycol		40 g/L
Sodium azide		0.95 g/L
R2	CRP TIA	Anti-CRP
Phosphate buffered saline		pH 7.43
Polyclonal goat anti-human CRP		(variable)
Sodium azide		0.95 g/L

These reagents are not classified as dangerous according to 1272/2008/EC regulation.

## REAGENTS PREPARATION

Ready for use.

## MATERIAL REQUIRED BUT NOT PROVIDED

1. Medical analysis laboratory equipment.
2. Spectrophotometer or Biochemistry Clinical Analyzer

## SAFETY CAUTIONS

- BIOLABO reagents are designated for professional use in laboratory
- Refer to current Material Safety Data Sheet available on request or on www.biolabo.fr
- Verify the integrity of the contents before use.
- Waste disposal: Respect legislation in force in the country.
- All specimens or reagents of biological origin should be handled as potentially infectious. Respect legislation in force in the country.
- ! Any serious incident that has occurred in connection with the device is notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is based.

## STABILITY AND STORAGE

Stored away from light, well cap in the original vial at 2-8°C, reagent is stable when stored and used as described in the insert:

- Unopened,
  - Until the expiry date stated on the label of the Kit.
- Once opened:
  - When free from contamination, 2 separated reagents are stable for :
    - 3 months at 2-8° C,
    - 24 h at room temperature
    - 30 days on the board of analysers.

## SPECIMEN COLLECTION AND HANDLING (5)

Use fresh serum.

If the test cannot be carried out on the same day, the serum may be stored at 2-8°C for 48 hours.

If stored for a longer period, the sample should be frozen.

## LIMITATIONS (6)

For a more comprehensive review of factors affecting this assay refer to the publication of Young D.S.

## CALIBRATION

- REF CRP CALSET51: CRP Standard Set
- Or
- REF CRP CALSH1: CRP Standard Super High (successive 1:2 dilutions in saline up to 6 different levels including zero point to generate calibration curve).
- Use saline as zero point

Calibration values are traceable to a reference material (RPPHS/CRM470) from the International Federation of Clinical Chemistry (IFCC).

The calibration frequency depends on proper instrument functions and on the preservation of the reagent.

Make a new calibration when changing reagent batch, if quality control results are found out of the confidence interval and after maintenance operations

## QUALITY CONTROL

- **REF** CRP CONTL1, **REF** CRP CONTL5: CRP Control Low
- **REF** CRP CONTH1, **REF** CRP CONTH5: CRP Control High
- External quality control program.

It is recommended to control in the following cases:

- At least once a run.
- At least once within 24 hours.
- When changing vial of reagent.
- After maintenance operations on the instrument.

If control is out of range, apply following actions:

1. Prepare a fresh control serum and repeat the test.
2. If control is still out of range, use a new vial of calibrator or a fresh calibrator and repeat the test.
3. If control is still out of range, repeat the tests with a new vial of reagent.

If control is still out of range, please contact BIOLABO technical support or your local Agent.

## REFERENCE INTERVAL (1)

IFCC Value: < 0.5 mg/dL

These values are applicable only to adults between 20 and 60 years of age

Each laboratory should establish its own normal ranges for the population that it serves.

## PERFORMANCES

On a clinical chemistry analyzer (Cobas Mira).

Detection limit: approximately 0.5 mg/dL

Linearity range: between 0.5 to 22 mg/dL.

Above 22 mg/dL, dilute the specimen with saline and re-assay taking into account the dilution factor to calculate the result. Linearity limit depends on specimen/reagent ratio.

Precision:	Low	Medium	High
<b>[%CV]</b>			
<b>Intra-Run</b>	4.06	2.57	3.44
<b>Inter-Run</b>		4.29	6.60

Accuracy:	Assigned	Measured
<b>[mg/dL]</b>		
BIOLABO	1.92 (1.63-2.21)	2.09
Biorad 2	2.59 (2.07-3.10)	2.78

Sensitivity: 0.0094 Abs/concentration unit

Specificity: Monospecific

Hook effect: > 84 mg/dL

Comparison study with Nephelometry:

$$y = 0.9981x - 0.0142 \quad r = 0.9917$$

Interferences:

No interference for: Haemoglobin (250 mg/dL), Na-citrate (1000 mg/dL), Heparin (50 mg/dL), Bilirubin (20 mg/dL), Triglyceride (2500 mg/dL)

Other substances may interfere with the result (see §Limitations)

## PROCEDURE

Let stand reagents, standards, control and specimens at room temperature.

Before use, mix reagent R2 by gentle swirling.

Manual Procedure:

Realise standard curve (§ Calibration)

Pipette into well identified test tubes:	Blank	Standards	Assays
<b>Buffer (R1)</b>	1000 µL	1000 µL	1000 µL
<b>Saline</b>	60 µL		
<b>Standards</b>		60 µL	
<b>Specimen</b>			60 µL

Mix well. Record absorbance A1 against blank at 340 nm

<b>Anti-CRP (R2)</b>	Blank	Standards	Assays
	100 µL	100 µL	100 µL

Mix and let stand for 5 minutes at room temperature.  
Record absorbance A2 against blank at 340 nm.

1- With Manual Procedure on Spectrophotometer, performances and stability data should be validated by user

2- Applications proposal are available on request of other analysers

## CALCULATION

Manual procedure:

Calculate  $\Delta Abs$  (Abs A2 – Abs A1) for standards, controls and assays.

Plot a Standard Curve "Concentration = f( $\Delta Abs$ )".

Read the concentration of controls and samples on the graph.













Automatic Biochemistry analyzer:

The analyzer provides directly final result.

For more details about calibration and calculation of results, refer to User's manual and specific application.

## REFERENCES

- (1) TIETZ N.W. *Text book of clinical chemistry*, 3<sup>rd</sup> Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p.493, p.481.
- (2) Marrack, J.R. and Richards, CB., *J. Immunol.* **20**, 1019 – 1040 (1971)
- (3) Ritchie, RF., *J. Lab. Clin. Med.* **70**, 512 - 517(1967)
- (4) Pepys MB. et al., *Ann. NY Acad. Sci.* **389**, 459 (1982)
- (5) *Clinical Guide to Laboratory Test*, 3<sup>rd</sup> Ed., N.W. TIETZ (1995) p. 919
- (6) YOUNG D.S., *Effect of Drugs on Clinical laboratory Tests*, 4<sup>th</sup> Ed. (1995) p. 3-511 to 3-512

 Manufacturer	 Expiry date	 In vitro diagnostic	 Storage temperature	 Dematerialized water	 Biological risk
 Product Reference	 See Insert	 Batch number	 Store away from light	 Sufficient for	 Dilute with