



BIOLABO
www.biolabo.fr

MANUFACTURER:
BIOLABO SAS,
Les Hautes Rives
02160, Maizy, France

Rheumatoid Factor Turbidimetric Immunoassay

Reagent for quantitative determination
of Rheumatoid Factor (RF) in human serum

REF RF050E	R1 1 x 50 mL	R2 1 x 10 mL
REF RF520E	R1 5 x 20 mL	R2 1 x 20 mL



Made in France

I: corresponds to significant modifications

TECHNICAL SUPPORT AND ORDERS

Tel: (33) 03 23 25 15 50

Fax: (33) 03 23 256 256

support@biolabo.fr

Latest revision: www.biolabo.fr

INTENDED USE

This quantitative test is intended to detect and measure rheumatoid factors in human serum. This test may be used with manual procedure on spectrophotometer or with Biochemistry Clinical Analyzer.

SUMMARY AND EXPLANATIONS (2) (4) (5)

RF is found in Rheumatoid Arthritis (RA), Sjogren's syndrome, scleroderma, sarcoidosis and other pathologies.

RF consists in a variety of antibodies that are usually of IgM class and will react with modified human IgG or IgG of animal origin.

About 75% of patients with RA have detectable RF of IgM class. The highest titres are seen in severe active chronic disease with vasculitis, and subcutaneous nodules.

80-90% of patients with Sjogren's syndrome will have high titre of RF.

PRINCIPLE (3)

Turbidimetric Immunoassay (TIA): Photometric measurement of antigen-antibody reaction between aggregated human IgG and RF, by the end-point method at 340 nm.

REAGENTS

R1	RF TIA	Buffer
Good's buffer (pH 7,4)		50 mM
Sodium azide		0.95 g/L
R2	RF TIA	Anti-RF
Heat-aggregated human IgG		(≤ 0,5 mg/mL)
Sodium azide		0.95 g/L

According to 1272/2008 regulation, these reagents are not classified as dangerous.

REAGENTS PREPARATION

Ready for use.

MATERIAL REQUIRED BUT NOT PROVIDED

1. Medical analysis laboratory equipment.
2. Spectrophotometer or Biochemistry Clinical Analyzer

SAFETY CAUTIONS

- BIOLABO reagents are designated for professional use in laboratory
- Refer to current Material Safety Data Sheet available on request or on www.biolabo.fr
- Verify the integrity of the contents before use.
- Waste disposal: Respect legislation in force in the country.
- All specimens or reagents of biological origin should be handled as potentially infectious. Respect legislation in force in the country.

I Any serious incident that has occurred in connection with the device is notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is based.

STABILITY AND STORAGE

Stored away from light, well cap in the original vial at 2-8°C, reagent is stable when stored and used as described in the insert:

Unopened,

- Until the expiry date stated on the label of the Kit.

Once opened:

- When free from contamination, 2 separated reagents are stable for:
 - 3 months at 2-8° C
 - 30 days on the board of analysers

SPECIMEN COLLECTION AND HANDLING (1) (2)

Use fresh serum.

If the test cannot be carried out on the same day, the serum may be stored at 2-8°C for 24 hours.

If stored for a longer period, the sample should be frozen once only.

LIMITATIONS (6)

For a more comprehensive review of factors affecting this assay refer to the publication of Young D.S.

CALIBRATION

- **REF** RF CALSET51: RF Standard Set
- Or
- **REF** RF CALSH1: RF Standard Super High (successive 1:2 dilutions in saline up to 5 different levels including zero point to generate calibration curve).
- Use saline as zero point

Calibration values are based on WHO Standardisation.

The calibration frequency depends on proper instrument functions and on the preservation of the reagent.

Make a new calibration when changing reagent batch, if quality control results are found out of the confidence interval and after maintenance operations

QUALITY CONTROL

- **REF** RF CONT1, **REF** RF CONT5: RF Control
- External quality control program.

It is recommended to control in the following cases:

- At least once a run.
- At least once within 24 hours.
- When changing vial of reagent.
- After maintenance operations on the instrument.

If control is out of range, apply following actions:

1. Prepare a fresh control serum and repeat the test.
2. If control is still out of range, use a new vial of calibrator or a fresh calibrator and repeat the test.
3. If control is still out of range, repeat the tests with a new vial of reagent.

If control is still out of range, please contact BIOLABO technical support or your local Agent.

REFERENCE INTERVAL (1)

WHO Values: 0-20 IU/mL

This range is given as a guide only.

Each laboratory should establish its own normal ranges for the population that it serves.

PERFORMANCES

On a clinical chemistry analyzer (Cobas Mira and Selectra 2).

Detection limit: 3 IU/mL

Linearity range: from 3 to 500 IU/mL.

Above 500 IU/mL, dilute the specimen with saline and re-assay taking into account the dilution factor to calculate the result. Linearity limit depends on specimen/reagent ratio.

	Low	Medium	High
Precision: [Intra-Run %CV]	2.68	1.38	1.55
[Inter-Run %CV]	3.07	1.40	1.78

	Assigned	Measured
Accuracy: [BIOLABO IU/mL]	111 (94-128)	105
Biorad 1	19.6 (16.6-22.5)	18.0
Biorad 2	39.8 (33.8-46.2)	39.7

Sensitivity: 0.00027 Abs/concentration unit

Specificity: Monospecific

Hook effect: No Risk

Comparison study with Nephelometry:

$$y = 0.6026 + 32. \quad r = 0.8776$$

Interferences:

No interference for: Haemoglobin (500 mg/dL), Ascorbic acid (50 mg/dL), Bilirubin (50 mg/dL), Intralipid (3%)

Other substances may interfere (see § Limitations)

PROCEDURE

Let stand reagents, standards, control and specimens at room temperature.

Before use, mix reagent R2 by gentle swirling.

Manual Procedure:

Realise calibration curve (§ Calibration)

Pipette into well identified test tubes:	Blank	Standards	Assays
Buffer (R1)	900 µL	900 µL	900 µL
Saline	50 µL		
Standards		50 µL	
Specimen			50 µL

Mix well. Record absorbance A1 against Blank at 340 nm.

Anti-RF (R2)	180 µL	180 µL	180 µL
---------------------	--------	--------	--------

Mix and let stand for 5 minutes at room temperature.
Record absorbance A2 against Blank at 340 nm.

1- With Manual Procedure on Spectrophotometer, performances and stability data should be validated by user

2- Applications proposal are available on request of other analyzers

CALCULATION

Manual procedure:

Calculate ΔAbs (Abs A2 – Abs A1) for standards, controls and assays. Plot a Standard Curve "Concentration = f(ΔAbs)".

Read the concentration (IU/mL) of controls and samples on the graph.

Automatic Biochemistry analyzer:

The analyzer provides directly final result.

For more details about calibration and calculation of results, refer to User's manual and specific application.

REFERENCES

- (1) TIETZ N.W. *Textbook of clinical chemistry*, 3rd Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p.1833.
- (2) *Clinical Guide to Laboratory Test*, 4th Ed., N.W. TIETZ (2006) p.958-961.
- (3) Klauss K., Bandilla M. D., and Mc Duffie M. D., *Arthritis and Rheumatism*, vol.12, n°2, p.74-81 (1969)
- (4) Waaler E., *Acta Path. Microb. Scan.*, 17 (1940)
- (5) Müller W., *The serology of Rheumatoid Arthritis*. Berlin-Göttingen-Heidelberg 97 (1962)
- (6) YOUNG D.S., *Effect of Drugs on Clinical laboratory Tests*, 4th Ed. (1995) p. 3-511 to 3-512

 Manufacturer	 Expiry date	 In vitro diagnostic	 Storage temperature	 Dematerialized water	 Biological risk
 Product Reference	 See Insert	 Batch number	 Store away from light	 Sufficient for	 Dilute with