



BIOLABO
www.biolabo.fr
MANUFACTURER:
BIOLABO SAS,
Les Hautes Rives
02160, Maizy, France

TRIGLYCERIDES GPO Method

Reagent for quantitative determination of triglycerides
in human serum and plasma

REF 80019	R1 2 x 50 mL	R2 2 x 50 mL	R3 1 x 5 mL
REF 87319	R1 10 x 100 mL	R2 10 x 100 mL	R3 1 x 5 mL



Made In France

TECHNICAL SUPPORT AND ORDERS

Tel : (33) 03 23 25 15 50

support@biolabo.fr

Latest revision : www.biolabo.fr

I: corresponds to significant modifications

INTENDED USE (1)


This reagent is designated for professional use in laboratory (manual or automated method).
It allows the quantification of triglycerides in human serum and plasma to complete lipids profile.

I GENERALITIES

The increase of triglycerides in blood can be of genetic origin or secondary to other metabolic disorders such as: diabetes mellitus, hyper and hypothyroidisms, hepatic diseases, acute and chronic pancreatitis, nephrosis. A rise in triglycerides also represents an atherogenic risk factor. It is responsible for the opalescence, or even the cloudiness of the serum. Corticoids and oestrogen/progestin treatments can also aggravate hypertriglyceridemia.

PRINCIPLE (4) (5)

Fossati and Prencipe method associated with Trinder reaction.
Reaction scheme is as follows:



The absorbance measured at 500 (480 – 520) nm is proportional to the amount of triglycerides in the specimen.

REAGENTS

R1 TRIGLYCERIDES GPO	Buffer
PIPES	100 mmol/L
Magnesium chloride	9.8 mmol/L
Chloro-4-phenol	3.5 mmol/L
Preservative	

According to 1272/2008 regulation, this reagent is not classified as dangerous

R2 TRIGLYCERIDES GPO	Enzymes
Lipase	≥ 1000 IU/L
Peroxidase (POD)	≥ 1700 IU/L
Glycerol 3 phosphate oxidase (GPO)	≥ 3000 IU/L
Glycerol Kinase (GK)	≥ 660 IU/L
4 - Amino – antipyrine (PAP)	0.5 mmol/L
Adenosine triphosphate Na (ATP)	1.3 mmol/L

Before reconstitution:

ATTENTION: Acute tox.4: H302 - Harmful if swallowed

P264 : Wash hands thoroughly after handling.

P301+312 : IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. P330 : Rinse mouth, P501: Dispose of contents/container in accordance with dangerous goods regulations. Classification due to : 4-Amino-antipyrine 1 - <2.5%. For more details refer to current Material Safety Data Sheet (MSDS)

Once reconstituted: Working Reagent is not classified as dangerous

R3 TRIGLYCERIDES GPO Standard 200 mg/dL (2.28 mmol/L)

According to 1272/2008 regulation, this reagent is not classified as dangerous

SAFETY CAUTIONS

BIOLABO reagents are designated for professional, in vitro diagnostic use (do not pipette with mouth).

- Refer to current Material Safety Data Sheet available on request or on www.biolabo.fr
- Verify the integrity of the contents before use.
- Waste disposal: Respect legislation in force in the country.

All specimens or reagents of biological origin should be handled as potentially infectious. Respect legislation in force in the country.

REAGENTS PREPARATION

Use a non-sharp instrument to remove aluminium cap.

Working reagent: Add promptly the contents of vial R2 into vial R1.

Mix gently and wait for complete dissolution.

Vial R3: Ready for use.

STABILITY AND STORAGE

Stored away from light, well capped in the original vial at 2-8°C, when stored and used as described, reagents are stable:

Unopened:

- Until expiry date stated on the label of the kit.

Once opened:

- Reconstitute immediately contents of vial R2.

- Standard (vial R3) :

Transfer requested quantity and store the vial at 2-8°C.

Once reconstituted:

- Transfer requested quantity and store in the original vial at 2-8°C.
- Working reagent is stable for 1 year when free from contamination.
- Discard reagent if cloudy or if absorbance at 500 nm > 0.200.
- Don't use working reagent after expiry date.

SPECIMEN COLLECTION AND HANDLING (2)

Serum or plasma (Heparin or EDTA) fasting ≥ 12 hours.

Separate from cells within 2 hours.

Do not use oxalate, fluoride or citrate.

Triglycerides are stable in specimen for:

- 5-7 days at 2-8°C.
- 3 months at -20°C.
- many years at -70°C.

Avoid repeated freezing and thawing.

LIMITS (1) (2) (3)

For a more comprehensive review of factors affecting this assay refer to the publication of Young D.S.

MATERIAL REQUIRED BUT NOT PROVIDED

1. Basic medical analysis laboratory equipment.
2. Spectrophotometer or Biochemistry Clinical Analyzer

QUALITY CONTROL

- **REF** 95010 EXATROL-N Level I
- **REF** 95011 EXATROL-P Level II
- External quality control program.

It is recommended to control in the following cases:

- At least once a run.
- At least once within 24 hours.
- When changing vial of reagent.
- After maintenance operations on the instrument.

If control is out of range, apply following actions:

1. Prepare a fresh control and repeat the test.
2. If control is still out of range, use a new vial of calibrator or a fresh calibrator and repeat the test.
3. If control is still out of range, repeat the tests with a new vial of reagent.

If control is still out of range, please contact BIOLABO technical support or your local Agent.

REFERENCE INTERVAL (6)

Serum or Plasma	mg/dL	[mmol/L]
Reference range	35-160	[0.40-1.82]

Each laboratory should establish its own normal ranges for the population that it serves.

PERFORMANCES

On Kenza 240TX, 37°C, 505nm

Detection limit: approx. 1 mg/dL

Precision:

Within-run N = 20	Low level	Normal level	High level	Between run N = 20	Low level	Normal level	High level
Mean mg/dL	41	102	206	Mean mg/dL	59	123	223
S.D. mg/dL	1.0	3.0	0.5	S.D. mg/dL	0.31	5.2	8.4
C.V. %	2.5	2.7	2.5	C.V. %	5.4	4.3	3.8

Analytical Sensitivity: approx. 0.125 abs at 500 nm for 100 mg/dL (manual procedure, 500 m, 1cm pathlength)

On Cobas Mira, 505 nm, 37°C:

Measuring Range: between 10 and 700 mg/dL

Comparison studies with commercially available reagent:

Realised with human specimens (n=75) between 25 and 350 mg/dL

$y = 1.0182x - 0.0302$ $r = 0.9958$

Interferences:

Ascorbic acid	Negative interference from 1 mg/dL
Total bilirubin	Positive interference from 100 µmol/L
Haemoglobin	No interference up to 300 µmol/L
Glucose	No interference up to 1260 mg/dL

Other substances may interfere (see § Limits)

CALIBRATION (7)

- **REF** 95015 Multicalibrator traceable to *SRM 909b2*

- Standard (vial R3): With manual procedure.

The calibration frequency depends on proper instrument functions and on the preservation of reagent.

Make a new calibration when changing reagent batch, if quality control results are found out of the range and after maintenance operations

PROCEDURE

Manual method

Let stand reagent and specimens at room temperature.

Reagent	1000 µL
Blank, Calibrator, control or specimen	10 µL
Mix well. Let stand for 10 minutes at room temperature or 5 minutes at 37°C. Record absorbance at 500 nm (480-520) against reagent blank. Color is stable for 1 h.	

- 1- Performances with manual procedure should be validated by user.
- 2- Kenza applications and other applications proposal are available on request.













CALCULATION

Serum or plasma:

$$\text{Result} = \frac{\text{Abs (Assay)}}{\text{Abs (Calibrator)}} \times \text{Calibrator concentration}$$

REFERENCES

- (1) *TIETZ N.W. Text book of clinical chemistry, 3rd Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p. 809-857.*
- (2) *Clinical Guide to Laboratory Test, 4th Ed., N.W. TIETZ (2006) p. 1074-1077.*
- (3) *YOUNG D.S., Effect of Drugs on Clinical laboratory Tests, 4th Ed. (1995) p.3-573 to 3-589*
- (4) *Fossati P., Prencipe L., Clin. Chem. (1982), 28, p.2077-2080.*
- (5) *Trinder P. Ann. Clin. Biochem. (1969), 6, p.27-29.*
- (6) *TIETZ N.W. Text book of clinical chemistry, 2nd Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1994)p. 1030-1058 et p. 1073-1080.*
- (7) *SRM: Standard Reference Material ®*

 Manufacturer	 Expiry date	 In vitro diagnostic	 Storage temperature	 Dematerialized water	 Biological risk
 Product Reference	 See Insert	 Batch number	 Store away from light	 Sufficient for	 Dilute with