



BIOLABO
www.biolabo.fr

MANUFACTURER :
BIOLABO SAS,
Les Hautes Rives
02160, Maizy, France

POTASSIUM Enzymatic method

Reagent for quantitative determination of potassium ions
in human serum or plasma.

REF 90084 R1 1x120 mL R2 1x30 mL R3 1x3 mL R4 1x3 mL



Made In France

TECHNICAL SUPPORT AND ORDERS

Tel: (33) 03 23 25 15 50

support@biolabo.fr

Latest revision: www.biolabo.fr

I: corresponds to significant modifications

INTENDED USE

This reagent is designated for professional use in laboratory (manual or automated method).

It allows the quantitative determination of potassium ions in human serum and plasma to determine if there is a problem with electrolyte balance. As part of a routine health check-up, results may be used as a screening test, in conjunction with other clinical signs and laboratory data.

GENERALITIES (1)

Disturbances in concentration of potassium may be encountered in case of hypokalaemia (metabolic alkalosis, metabolic acidosis, and abnormal acid-base balance), hyperkalaemia (over-administration of potassium, acidosis or crush injuries), kidneys damage, Addison disease or other diseases involved in electrolytes imbalance.

PRINCIPLE (1) (2)

Potassium is determined spectrophotometrically through a kinetic coupling assay system using potassium dependent pyruvate kinase. Pyruvate generated is converted to lactate accompanying conversion of NADH in $\text{NAD}^+ + \text{H}^+$. The corresponding decrease of optical density at 380 nm is proportional to the potassium concentration in the serum.

REAGENTS

R1 Potassium	Reagent 1		
LDH	< 50	KU/L	
NADH analog substrate	< 10	mmol/L	
Sodium Azide	0.05	%	
R2 Potassium	Reagent 2		
Pyruvate Kinase	< 50	KU/L	
Sodium Azide	0.05	%	
R3 Potassium	Cal 1	Calibrator	
Level 1 Potassium Chloride	approx. 3	mmol/L	
Sodium azide	< 0.1	%	
R4 Potassium	Cal 2	Calibrator	
Level 2 Potassium Chloride	approx. 7	mmol/L	
Sodium azide	< 0.1	%	

According to 1272/2008/EC Regulation, these reagents are not classified as dangerous.

SAFETY CAUTIONS

- Refer to current Material Safety Data Sheet available on request or on www.biolabo.fr
 - Verify the integrity of the contents before use.
 - Waste disposal: Respect legislation in force in the country.
 - All specimens or reagents of biological origin should be handled as potentially infectious. Respect legislation in force in the country.
- Any serious incident that has occurred in connection with the device is notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is based.

REAGENTS PREPARATION

Ready for use.

STABILITY AND STORAGE

Stored away from light, well capped in the original vial at 2-8°C, when stored and used as described, reagents are stable:

Unopened:

- Until expiry date stated on the label of the kit.

Once opened:

- Reagents are stable for at least 30 days.

SPECIMEN COLLECTION AND HANDLING (3)(4)

Plasma (Lithium heparinate)

- Do not use plasma from blood that has been stored in ice water
- Remove from cells promptly and test as soon as possible after collection.

Unhaemolyzed serum:

- Separate from cells as soon as possible after collection.

LIMITS (5)

For a more comprehensive review of factors affecting this assay refer to the publication of Young D.S.

MATERIALS REQUIRED BUT NOT PROVIDED

- Basic medical analysis laboratory equipment.
- Spectrophotometer or Biochemistry Clinical Analyzer

Manufacturer	Expiry date	In vitro diagnostic	Storage temperature	Dematerialized water	Biological risk
Product Reference	See Insert	Batch number	Store away from light	Sufficient for	Dilute with

QUALITY CONTROL

- **REF** 95010 EXATROL-N Level I
- **REF** 95011 EXATROL-P Level II
- External quality control program

It is recommended to control in the following cases:

- At least once a run
- At least once within 24 hours
- When changing vial of reagent
- After maintenance operations on the instrument

If control is out of range, apply following actions:

1. Prepare a fresh control serum and repeat the test
 2. If control is still out of range, use a new vial of fresh calibrator
 3. If control is still out of range, use a new vial of reagent and reassay
- If control is still out of range, please contact BIOLABO technical support or your local Agent.

REFERENCE INTERVAL (3) (4)

Serum or plasma	mEq/L	[mmol/L]
Premature, cord	5.0-10.2	[5.0-10.2]
Premature, 48h	3.0-6.0	[3.0-6.0]
New born, in cord	5.6-12.0	[5.6-12.0]
New Born	3.7-5.9	[3.7-5.9]
Infant	4.1-5.3	[4.1-5.3]
Child	3.4-4.7	[3.4-4.7]
Thereafter	3.5-5.1	[3.5-5.1]
Men	3.5-4.5	[3.5-4.5]
Women	3.4-4.4	[3.4-4.4]

Each laboratory should establish its own normal ranges for the population that it serves.

PERFORMANCES

On automatic analyzer AU400, at 37°C, 380 nm

Linearity Range: between 2 mmol/L and 8 mmol/L

Precision:

Within-run N = 20	Low level	Normal level	Between run N = 20	Low level	Normal level
Mean (mmol/L)	4.62	6.96	Mean (mmol/L)	4.62	6.96
S.D. mmol/L	0.052	0.084	S.D. mmol/L	0.081	0.122
C.V. %	1.1	1.2	C.V. %	1.8	1.8

Analytical Sensitivity:

approx. 0.101 abs / 1mmol/L (380 nm, 1 cm path length, 37°C)

Interferences (less than 5% deviation for listed concentrations):

Interference	Concentration	Interference	Concentration
Ascorbic acid	10 mmol/L	NH4+	0,5 mmol/L
Triglycerides	1000 mg/dL	Ca2+	7,5 mmol/L
Hemoglobin	500 mg/dL	Pi (inorganic phosphorous)	2 mmol/L
Conjugate bilirubin	20 mg/dL	Fe3+	0,5 mmol/L
Bilirubin	15 mg/dL	Cu2+	0,5 mmol/L
Na+	150 mmol/L	Zn2+	0,5 mmol/L

Other substances may interfere (see § Limits)

Comparison studies with commercially available reagent:

Automated analyzer (specimens n=56) from 2.5 to 7.8 mmol/L

$$y = 1,07 x - 0,3042$$

$$r = 0,9902$$

On the board stability: 2 separate reagents are stable 30 days.

Calibration Frequency: 14 days

Make a new calibration when changing reagent batch, if quality control results are found out of the established range and after maintenance operations.

CALIBRATION (6)

- Cal1 and Cal 2 (Vial R3 and R4) traceable to SRM956

The calibration frequency depends on proper instrument functions and on the preservation of reagent.

PROCEDURE

Manual method

Let stand reagent and specimens at room temperature.

Reagent 1	800 µL
Blank, Standards, control or specimen	20 µL
Mix well. Let stand for 5 minutes at 37°C	
Reagent 2	200 µL
Mix well. Read at 380 nm absorbance A1 after 60 sec and A2 after 240 sec . Calculate ΔAbs (Abs A2 – Abs A1) for Blank, Standards and Assays.	

- 1- Performances with manual procedure should be validated by user.
- 2- KENZA applications and other applications proposal are available on request.

CALCULATION

Serum or plasma:

$$\frac{\Delta \text{Abs (Assay)} - \Delta \text{Abs (Blank)}}{\Delta \text{Abs (Standard)} - \Delta \text{Abs (Blank)}}$$

$$\Delta \text{Abs (Standard)} - \Delta \text{Abs (Blank)}$$

Interpolate the ΔA in the Calibration Curve

REFERENCES

- (1) Bergmeyer, H.U., Gawehn, K., and Grassl, M. (1974) in *Methods of Enzymatic Analysis. Second Edition, Volume I*, 509-510, Academic Press, Inc., New York.
- (2) M.N. Berry, R. D. Mazzachi, M. Pejakovic, and M. J. Peake *Enzymatic Determination of Potassium in Serum. CLIN. CHEM.* 35/5, 817-820 (1989).
- (3) Tietz N.W. *Text book of clinical chemistry*, 3rd Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p. 1058, 1101—1104.
- (4) *Clinical Guide to Laboratory Test*, 4th Ed., N.W. Tietz (2006) p. 880.
- (5) YOUNG D.S., *Effect of Drugs on Clinical laboratory Tests*, 4th Ed. (1995) p. 3-476 à 3-486
- (6) SRM: Standard Reference Material®