

BIOLABO www.biolabo.fr

MANUFACTURER: **BIOLABO S.A.S** Les Hautes Rives 02160, Maizy, France

	D1 4 x 17 ml
REFRIDUZ	K1 4X1711L
REF K2002	R1 4 x 25 mL

ALBUMIN BCG Method

Reagent for quantitative determination of albumin in human serum and plasma.

TECHNICAL SUPPORT AND ORDERS

((

Made In France



I INTENDED USE

Tel : (33) 03 23 25 15 50

Latest revision : www.biolabo.fr

support@biolabo.fr

Reagent for quantitative determination of albumin in human serum or plasma to screen its level.

For professional laboratory Use (automated method).

GENERALITIES (3)

Albumin is the most abundant plasma protein. The primary function of albumin is generally considered to be the maintenance of colloidal osmotic pressure (COP) in both the vascular and extravascular spaces. Albumin have the ability to bind and transport a large number of compounds such as free fatty acids, phospholipids, metallic ions, amino acids, drugs, hormones, bilirubin, among many others.

PRINCIPLE (1) (2)

In buffered solution at pH 4.2, bromocresol green binds albumin to form a coloured compound which absorbance, measured at 630 nm (620-640) is proportional to the albumin concentration in the specimen.

REAGENT COMPOSITION

R1	ALB	Reagent

Succinic acid	83	mmol/l	_
Bromocresol green (BCG)		167	µmol/L
Sodium hydroxide		50	mmol/L
Polyoxyethylene monolauryl e	ther	1.00	g/L
Preservative			

ATTENTION, Met. Corr.1: H290 - May be corrosive for metals

P234: Keep only in original container, P390: Absorb spillage to prevent material damage. Classification due to: Sodium Hydroxide < 1% For more details, refer to Safety data sheet (SDS)

SAFETY CAUTIONS

- · Refer to current Material Safety Data Sheet available on request or on www.biolabo.fr
- · Verify the integrity of the contents before use.
- Waste disposal: Respect legislation in force in the country.
- All specimens or reagents of biological origin should be handled as potentially infectious. Respect legislation in force in the country.

Any serious incident that has occurred in connection with the device is notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is based.

REAGENT PREPARATION

Ready for use.

STABILITY AND STORAGE

Stored away from light, well cap in the original vial at 2-8°C, reagent is stable when stored and used as described in the insert:

Unopened,

• Until the expiry date stated on the label of the Kit.

Once opened:

- · Reagents are stable at least 1 year.
- Discard reagent if cloudy or if absorbance at 630 nm is > 0.300.

SPECIMEN COLLECTION AND HANDLING (4)

Serum or plasma. Serum must be separated from blood cells within 2 hours.

- Serum albumin is stable in serum for:
 - √72 hours at 2-8°C
 - ✓ 6 months at 20°C.

LIMITS (4) (5) (6) (7)

Heparinised plasma gives higher values than serum.

This interference can be avoided by working with dichromatic procedure (2nd wavelength is 550 nm or 700 nm).

Clofibrate and Phenylbutazone decrease albumin value with this procedure.

For a more comprehensive review of factors affecting this assay refer to the publication of Young D.S.

MATERIAL REQUIRED BUT NOT PROVIDED

- 1.Basic medical analysis laboratory equipment.
- 2. Biochemistry Clinical Analyzer Kenza One, Kenza 240TX/ISE or Kenza 450TX/ISE

REFERENCE INTERVAL (4)

Albumin	g/dL	[µmol/L]
0 to 4 days	2.8-4.4	[421-662]
4 days to 14 years	3.8-5.4	[572-813]
14 to 18 years	3.2-4.5	[482-677]
18 to 60 years	3.4-4.8	[512-722]
60 to 90 years	3.2-4.6	[482-692]
> 90 years	2.9-4.5	[436-677]

Each laboratory should establish its own normal ranges for the population that it serves.

PERFORMANCES at 37°C on KENZA ONE

Linearity Range: between 0.53 g/dL and 7.0 g/dL Detection limit: approx. 0.01 g/dL

Precision:

Within-run N = 20	Low level	Normal level	High Ievel	Between run N = 20	Low level	Normal level	High level
Mean (g/dL)	2.15	3.11	4.61	Mean (g/dL)	2.12	2.99	4.55
S.D. g/dL	0.042	0.065	0.046	S.D. g/dL	0.035	0.044	0.068
C.V. %	2.0	2.1	1.0	C.V. %	1.7	1.5	1.5

Comparison studies with commercially available reagent: In clinical environment with specimens between 1.1 and 4.5 g/dL (n=108)

y = 0.9028 x + 0.2441R= 0.9865

Analytical Sensitivity: approx. 0.0125 abs (0.1g/dL)

Interferences:

Turbidity	Positive interference from 0.295 OD
Total bilirubin	No interference up to 533 µmol/L
Direct bilirubin	No interference up to 486 µmol/L
Ascorbic acid	No interference up to 2500 mg/dL
Glucose	No interference up to 990 mg/dL
Haemoglobin	Positive interference from 333 µmol/L

Other substances may interfere (see § Limits)

I On the board stability: 2 months

Calibration Stability: 1 month

Make a new calibration when changing reagent batch, if quality control results are found out of the established range and after maintenance operations.

Performances and stability data on Kenza 450TX/ISE and Kenza 240TX/ISE are available on request

CALIBRATION (8)

REF BIOLABO Multicalibrator traceable to SRM 927

The calibration frequency depends on proper instrument functions and on preservation of the reagent.

QUALITY CONTROL

- REF 95010 EXATROL-N Level I
- REF 95011 EXATROL-P Level II
- External quality control program.
 - It is recommended to control in the following cases:
 - At least once a run
 - At least once within 24 hours
 - When changing vial of reagent
 - · After maintenance operations on the instrument
 - If control is out of range, apply following actions:
 - 1. Prepare a fresh control serum and repeat the test
 - 2. If control is still out of range, use a new vial of fresh calibrator

3. If control is still out of range, use a new vial of reagent and reassay If control is still out of range, please contact BIOLABO technical support or your local Agent

PROCEDURE

Refer to validated application of the Kenza Analyzer used

CALCULATION

The analyzer provides directly final result.

Refer to the instruction of use of Kenza analyzer

REFERENCES

- Albumin standards and the measurement of serum albumin with bromcresol (1) green, DOUMAS B.T., WATSON W.A., BIGGS H.G. - Clin. Chim. Acta., 31, (1971), p. 87-96.
- Determination of serum albumin, DOUMAS B.T. and BIGGS H.G. -(2)Standard methods of clinical chemistry - Acad. Press. N.Y. Vol 7 (1972) p. 175-188.
- (3) TIETZ N.W. Text book of clinical chemistry, 3rd Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p. 482-485. Clinical Guide to Laboratory Test, 4th Ed., N.W. TIETZ (2006) p. 68-71
- YOUNG D.S., Effect of Drugs on Clinical laboratory Tests, 4th Ed. (1995) (5) p.3-16 to 3-22
- Overestimation of Albumin in Heparinized Plasma, HALLBACH J., (6) HOFFMANN G.E., GUDER W.G., Clin. Chem. Vol 37 No 4 (1991), p. 566-568.
- Improved specificity of serum Albumin determination and estimation of (7)"acute phase reactants" by use of the bromcresol green reaction. Jan E. C. Gustafsson, Clin. Chem., Vol 22,n°5, (1976) p.616-622
- (8) SRM: Standard Reference Material ®

	Ω	IVD	X	H ₂ O	¢
Manufacturer	Expiry date	In vitro diagnostic	Storage temperature	Dematerialized water	Biological risk
REF		LOT	淡	E	\rightarrow
Product Reference	See Insert	Batch number	Store away from light	Sufficient for	Dilute with