



**BIOLABO**  
www.biolabo.fr  
**MANUFACTURER:**  
**BIOLABO SAS,**  
Les Hautes Rives  
02160, Maizy, France

I REF	K1002	R1	4 x 17 mL
I REF	K2002	R1	4 x 25 mL

# ALBUMIN BCG Method

Reagent for quantitative determination of albumin  
in human serum and plasma.

## TECHNICAL SUPPORT AND ORDERS

Tel : (33) 03 23 25 15 50  
support@biolabo.fr  
Latest revision : www.biolabo.fr



Made In France

I: corresponds to significant modifications



## INTENDED USE

This reagent is designated for professional use in laboratory (automated method). It allows the quantification of global activity of the alkaline phosphatase enzyme in human serum or plasma.

## GENERALITIES (3)

Albumin is the most abundant plasma protein. The primary function of albumin is generally considered to be the maintenance of colloidal osmotic pressure (COP) in both the vascular and extravascular spaces. Albumin have the ability to bind and transport a large number of compounds such as free fatty acids, phospholipids, metallic ions, amino acids, drugs, hormones, bilirubin, among many others.

## PRINCIPLE (1) (2)

In buffered solution at pH 4.2, bromocresol green binds albumin to form a coloured compound which absorbance, measured at 630 nm (620-640) is proportional to the albumin concentration in the specimen.

## REAGENT COMPOSITION

R1	ALB	Reagent	
		Succinic acid	83 mmol/L
		Bromocresol green (BCG)	167 µmol/L
		Sodium hydroxide	50 mmol/L
		Polyoxyethylene monolauryl ether	1.00 g/L
		Preservative	

**ATTENTION**, Met. Corr.1: H290 – May be corrosive for metals

P234: Keep only in original container, P390: Absorb spillage to prevent material damage. Classification due to: Sodium Hydroxide < 1% For more details, refer to Safety data sheet (SDS)

## SAFETY CAUTIONS

- Refer to current Material Safety Data Sheet available on request or on www.biolabo.fr
- Verify the integrity of the contents before use.
- Waste disposal: Respect legislation in force in the country.
- All specimens or reagents of biological origin should be handled as potentially infectious. Respect legislation in force in the country.

Any serious incident that has occurred in connection with the device is notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is based.

## REAGENT PREPARATION

Ready for use.

## STABILITY AND STORAGE

**Stored away from light, well cap in the original vial at 2-8°C, reagent is stable when stored and used as described in the insert:**

Unopened,

- Until the expiry date stated on the label of the Kit.

Once opened:

- Reagents are stable at least 1 year.
- Discard reagent if cloudy or if absorbance at 630 nm is > 0.300.

## SPECIMEN COLLECTION AND HANDLING (4)

Serum or plasma. Serum must be separated from blood cells within 2 hours.

Serum albumin is stable in serum for:

- ✓ 72 hours at 2-8°C.
- ✓ 6 months at -20°C.

## LIMITS (4) (5) (6) (7)

Heparinised plasma gives higher values than serum. This interference can be avoided by working with dichromatic procedure (2<sup>nd</sup> wavelength is 550 nm or 700 nm).

Clofibrate and Phenylbutazone decrease albumin value with this procedure.

For a more comprehensive review of factors affecting this assay refer to the publication of Young D.S.

## MATERIAL REQUIRED BUT NOT PROVIDED

1. Basic medical analysis laboratory equipment.
2. Biochemistry Clinical Analyzer Kenza One, Kenza 240TX/ISE or Kenza 450TX/ISE

## REFERENCE INTERVAL (4)

Albumin	g/dL	[ $\mu$ mol/L]
0 to 4 days	2.8-4.4	[421-662]
4 days to 14 years	3.8-5.4	[572-813]
14 to 18 years	3.2-4.5	[482-677]
18 to 60 years	3.4-4.8	[512-722]
60 to 90 years	3.2-4.6	[482-692]
> 90 years	2.9-4.5	[436-677]

Each laboratory should establish its own normal ranges for the population that it serves.

## PERFORMANCES at 37°C on KENZA ONE

Linearity Range: between 0.53 g/dL and 7.0 g/dL

Detection limit: approx. 0.01 g/dL

Precision:

Within-run N = 20	Low level	Normal level	High level	Between run N = 20	Low level	Normal level	High level
Mean (g/dL)	2.15	3.11	4.61	Mean (g/dL)	2.12	2.99	4.55
S.D. g/dL	0.042	0.065	0.046	S.D. g/dL	0.035	0.044	0.068
C.V. %	2.0	2.1	1.0	C.V. %	1.7	1.5	1.5

Comparison studies with commercially available reagent:  
In clinical environment with specimens between 1.1 and 4.5 g/dL  
(n=108)

$$y = 0.9028x + 0.2441 \quad R = 0.9865$$

Analytical Sensitivity: approx. 0.0125 abs (0.1g/dL)

Interferences:

Turbidity	Positive interference from 0.295 OD
Total bilirubin	No interference up to 533 $\mu$ mol/L
Direct bilirubin	No interference up to 486 $\mu$ mol/L
Ascorbic acid	No interference up to 2500 mg/dL
Glucose	No interference up to 990 mg/dL
Haemoglobin	Positive interference from 333 $\mu$ mol/L

Other substances may interfere (see § Limits)

On the board stability: 3 months

Calibration Stability: 1 month

Make a new calibration when changing reagent batch, if quality control results are found out of the established range and after maintenance operations.

Performances and stability data on Kenza 450TX/ISE and Kenza 240TX/ISE are available on request

## CALIBRATION (8)

- **REF** BIOLABO Multicalibrator traceable to SRM 927  
The calibration frequency depends on proper instrument functions and on preservation of the reagent.

## QUALITY CONTROL

- **REF** 95010 EXATROL-N Level I
  - **REF** 95011 EXATROL-P Level II
  - External quality control program.
- It is recommended to control in the following cases:

- At least once a run
  - At least once within 24 hours
  - When changing vial of reagent
  - After maintenance operations on the instrument
- If control is out of range, apply following actions:
1. Prepare a fresh control serum and repeat the test
  2. If control is still out of range, use a new vial of fresh calibrator
  3. If control is still out of range, use a new vial of reagent and reassay
- If control is still out of range, please contact BIOLABO technical support or your local Agent

## PROCEDURE

Refer to validated application of the Kenza Analyzer used













## CALCULATION

The analyzer provides directly final result.

Refer to the instruction of use of Kenza analyzer

## REFERENCES

- (1) *Albumin standards and the measurement of serum albumin with bromocresol green*, DOUMAS B.T., WATSON W.A., BIGGS H.G. - *Clin. Chim. Acta.*, **31**, (1971), p. 87-96.
- (2) *Determination of serum albumin*, DOUMAS B.T. and BIGGS H.G. - *Standard methods of clinical chemistry* - Acad. Press. N.Y. Vol 7 (1972) p. 175-188.
- (3) *TIETZ N.W. Text book of clinical chemistry*, 3<sup>rd</sup> Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p. 482-485.
- (4) *Clinical Guide to Laboratory Test*, 4<sup>th</sup> Ed., N.W. TIETZ (2006) p. 68-71
- (5) *YOUNG D.S., Effect of Drugs on Clinical laboratory Tests*, 4<sup>th</sup> Ed. (1995) p.3-16 to 3-22
- (6) *Overestimation of Albumin in Heparinized Plasma*, HALLBACH J., HOFFMANN G.E., GUDER W.G., *Clin. Chem.* Vol 37 No 4 (1991), p. 566-568.
- (7) *Improved specificity of serum Albumin determination and estimation of "acute phase reactants" by use of the bromocresol green reaction*. Jan E. C. Gustafsson, *Clin. Chem.*, Vol 22, n°5, (1976) p.616-622
- (8) *SRM: Standard Reference Material* ®

 Manufacturer	 Expiry date	 In vitro diagnostic	 Storage temperature	 Dematerialized water	 Biological risk
 Product Reference	 See Insert	 Batch number	 Store away from light	 Sufficient for	 Dilute with