



**BIOLABO**  
**www.biolabo.fr**  
**MANUFACTURER:**  
**BIOLABO S.A.S**  
Les Hautes Rives  
02160, Maizy, France

# EXTRA Cleaning

## Degreaser Solution

REF CO4015 2 x 250 mL

### TECHNICAL SUPPORT AND ORDERS

Tel: (33) 03 23 25 15 50

support@biolabo.fr

Latest revision: www.biolabo.fr



Made in France

I: corresponds to significant modifications

### INTENDED USE

This reagent is designated for professional use in laboratory (automated method). It allows the cleaning of reaction cuvettes of BIOLABO KENZA analysers to maintain them in perfect conditions.

### GENERALITIES (1) (2)

Periodic cleaning of reading analytical system, hydraulic circuits and reaction cuvettes is part of operations of ordinary preventive maintenance which must be performed by the end of user in order to achieve the best performance during analytical routine.

### I PRINCIPLE

This solution is to be used with IPO Cleaning. Washing with the two solutions is highly recommended as it has a deproteinizing/degreasing function and helps to maintain the necessary transparency of reaction cuvettes. From the wash station, the cleaning water flows to a liquid waste container.

### REAGENTS

<b>R1</b>	<b>EXTRA Cleaning</b>	Degreaser Solution
Sulfuric acid	0,075%	

According to CLP Regulation (EC) 1272/2008, this reagent is not classified as dangerous:

### MATERIAL REQUIRED BUT NOT PROVIDED

- KENZA ONE, KENZA 240TX/ISE, KENZA 450TX/ISE analyser.
- REF CO4020 IPO Cleaning

### SAFETY CAUTIONS

- Refer to current Material Safety Data Sheet available on request or on www.biolabo.fr
  - Verify the integrity of the contents before use.
  - Waste disposal: Respect legislation in force in the country.
  - All specimens or reagents of biological origin should be handled as potentially infectious. Respect legislation in force in the country.
- Any serious incident that has occurred in connection with the device is notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is based

### REAGENTS PREPARATION

This solution is ready to use.

### STABILITY AND STORAGE

**Stored away from light, well cap in the original vial at 15-25°C, reagent is stable when stored and used as described in the insert:**

Unopened,

- Until the expiry date stated on the label of the Kit.

Once opened:

- when free from contamination, stored and used as described in the insert, the contents of the vial may be used for at least 3 months.
  - ✓ Do not refill container
  - ✓ Discard any cloudy reagent
  - ✓ Don't use after expiry date stated on the label

### PERFORMANCES

No cross contamination detected using this Cleaning Solution if used as indicated in the User's manual §Maintenance.

To be used as indicated in the User's manual of KENZA Analyser used.

### REFERENCES

- (1) Occupational Safety and Health Standards; Bloodborne pathogens (29CFR1910.1030) Federal Register July 1, (1998) ; 6, p.267-280
- (2) Directive du conseil de l'Europe (90/679/CEE) J. O. de la communauté européenne n°L374 du 31.12.1990, p.1-12

Manufacturer	Expiry date	In vitro diagnostic	Storage temperature	Dematerialized water	Biological risk
Product Reference	See Insert	Batch number	Store away from light	Sufficient for	Dilute with